

## Product datasheet for MC225949

### Grin1os (NM\_001242940) Mouse Untagged Clone

#### Product data:

|                           |                                                                                    |
|---------------------------|------------------------------------------------------------------------------------|
| Product Type:             | Expression Plasmids                                                                |
| Product Name:             | Grin1os (NM_001242940) Mouse Untagged Clone                                        |
| Tag:                      | Tag Free                                                                           |
| Symbol:                   | Grin1os                                                                            |
| Synonyms:                 | E130003G02Rik                                                                      |
| Mammalian Cell Selection: | Neomycin                                                                           |
| Vector:                   | pCMV6-Entry (PS100001)                                                             |
| E. coli Selection:        | Kanamycin (25 ug/mL)                                                               |
| Fully Sequenced ORF:      | >MC225949 representing NM_001242940<br>Red=Cloning site Blue=ORF Orange=Stop codon |

TTTTGTATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCCCGCATCGCC

ATGGAAGGCGCACCTGCCATGTTCTCAAAGGGCTGTCTTGCGCATGCCGATGCCAAAGCCGGAGCGGA  
 AGAACAGCTACCCGTGGTCACCAGATCGCACTTCTGTGAAGCCTCAAACCTCCAGCACAGCTGAGTCCCA  
 GATGAAGGCATGGAGCTTGCTGCCGGCAGAGTCAGGGTCAAAGGAGCCCCGTGTAGCCTGTTACCCCCAC  
 CCTGCTCACACTGCAGAAGACAGTCTTCACTTGTGTTTGTGGTGAGAGGGCTCCCTTTTCCTGTCCC  
 TTTACTGTCCCCCAGCCAGGGCGGCCCTTGGTCTGCCTGCTCACTTGTCCCGCACAGCCTGGATGGCC  
 TCAGCTGCACTCTCATAATTGTGCTTCTCCATGTGCCGGTACATGGTGCTCAACTCCACCTGCCTCCGGA  
 AGTAGATATCCAGAGCTCTGTTTTACAGTTGCATAGATGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

|                    |              |
|--------------------|--------------|
| Restriction Sites: | Sgfl-MluI    |
| ACCN:              | NM_001242940 |
| Insert Size:       | 462 bp       |

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).


[View online »](#)

|                               |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
|-------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>OTI Annotation:</b>        | Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.                                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Components:</b>            | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).                                                                                                                                                                                                                                                                                                                              |
| <b>Reconstitution Method:</b> | <ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol> |
| <b>Note:</b>                  | Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>RefSeq:</b>                | <u><a href="#">NM_001242940.1</a></u> , <u><a href="#">NP_001229869.1</a></u>                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>RefSeq Size:</b>           | 2894 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>RefSeq ORF:</b>            | 462 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Locus ID:</b>              | 320354                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Cytogenetics:</b>          | 2 A3                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |