

Product datasheet for MC225983

7530416G11Rik (NM_001256072) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	7530416G11Rik (NM_001256072) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	7530416G11Rik
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>MC225983 representing NM_001256072 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGCTCACAAGGGCACTGCCTGGATTGCTGTCTACAGAGGCTCGGCTCCCTCTGATGTTTTGAACCCAG
 GGAATCACTGGCAAGAATTATTTTCAGCTCCCAAGAAGGGGACTTGGCGACAGACACAGAAGGAGCTGTG
 CGGGCAGAGAGGGCCAGCAGGGTTGCCATGGGAACGCCTCCTGATTGAAAGGGATTCTATTGGGGTGCG
 AAGGACAGAAGGAGGCACCCTGAATCCTGTCTACCTACCCATTCCCTCCCTCAGCCTGGGACCAACTCC
 AAGGCAGGCAGCCTCAGGCTTGCCAGGCTTCTGGAGGGCTGGAGAAGGAATTGGAATTGACATTGAATT
 CCATGTATGCTGGAGCCTGGTCTCTATTCCAGGCAGACAGTAGCAACCAAGCAAGCAAATGGATGGAGA
 AGGAAATACGGGGCCTACACACAATGGAGATTTACTGAGTCTCCAAGAGGAACGGGATGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	SgfI-MluI
ACCN:	NM_001256072
Insert Size:	480 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).


[View online »](#)

OTI Annotation:	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_001256072.1</u> , <u>NP_001243001.1</u>
RefSeq Size:	1161 bp
RefSeq ORF:	480 bp
Locus ID:	328577
Cytogenetics:	15 E2