

Product datasheet for **MC226017**

Jdp2 (NM_001205052) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Jdp2 (NM_001205052) Mouse Untagged Clone
Tag: Tag Free
Symbol: Jdp2
Synonyms: Jundm2; Jundp2; TIF
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC226017 representing NM_001205052
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGATGCCTGGGCAGATCCCAGACCTTCAGTGACCGCAGGCTCTCTGCCAGGGCTCGGCCCTCACCG
GACTTCCCAGCTCTGCTCTGACCACAGAGGAGCTGAAATACGCTGACATCCGCAACATTGGGGCGATGAT
TGCGCCCTTGCACTTCTGGAGGTGAACTGGGCAAGAGGCCCAACCCGTGAAGAGTGAGCTAGACGAG
GAAGAAGAGCGAAGGAAAAGGCGCCGGGAAAAGAACAAGTGCCTGCAGCCAGATGCCGGAACAAGAAGA
AGGAACGCACAGAGTTTCTGCAGAGGGAGTCAGAGCGGCTGGAGCTCATGAACGCAGAGCTGAAGACGCA
GATAGAGGAGCTGAAGCTGGAGCGGCAACAGCTTATCCTGATGCTCAACCGCCACCGCCCACTGCATC
GTGCGCACAGATAGCGTCAGGACGCCCGAGTCCGAAGGCAACCCACTGCTGGAGCAGCTGGACAAGAAGT
GA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001205052
Insert Size: 492 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

OTI Annotation:	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001205052.1</u> , <u>NP_001191981.1</u>
RefSeq Size:	1567 bp
RefSeq ORF:	492 bp
Locus ID:	81703
UniProt ID:	<u>P97875</u>
Cytogenetics:	12 D2
Gene Summary:	<p>Component of the AP-1 transcription factor that represses transactivation mediated by the Jun family of proteins. Involved in a variety of transcriptional responses associated with AP-1, such as UV-induced apoptosis, cell differentiation, tumorigenesis and antitumogenesis. Can also function as a repressor by recruiting histone deacetylase 3/HDAC3 to the promoter region of JUN. May control transcription via direct regulation of the modification of histones and the assembly of chromatin (By similarity).[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (2) differs in the 5' UTR, compared to variant 1. Variants 1, 2 and 3 encode the same protein. Sequence Note: The RefSeq transcript and protein were derived from genomic sequence to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on alignments.</p>