

Product datasheet for **MC226047**

Ssxb6 (NM_001205108) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Ssxb6 (NM_001205108) Mouse Untagged Clone
Tag: Tag Free
Symbol: Ssxb6
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Fully Sequenced ORF: >MC226047 representing NM_001205108
Red=Cloning site **Blue**=ORF **Orange**=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGAACAGTGAGTTCCTGTGAAAAGATCCCCATGGAAGTGCTTTATGAACAAAGAATATCTGCAAGG
CATTCCAGGATATTTCTACATACTTCTCTGATGAAGAGTGGGTAAAGCTGACTCAATGGCAGAAAAGTGC
CTATGTGTACATGAAAAGAACTACATCAGAATGACTGACCTAGGTGTCACTGTGAACCAACCAGTTTTC
ATGCGTGGCAAGGAGCAGGCCAAGCAATCCCTGGTGAAGGCATTGAAATTCATGACAGTGATGATGAAT
GCTTTGAAGGATCTTTTGGTGTGACACCGATAAAACGAATGAAGCTGACATCAGTGACAATAAGTATTCA
TAATGTAGAAGGAAGCCATGACTCTGGAGAGAATGATTCAAACCTGGCAGAACTGGTGGCATTTCAGGTC
AATGTTTGGAGCCACCGCTGCGAGAGAGGAAGTATCGAGTGATATATTCAGAGATTAGTGATCCTGAAG
AGGAAGAAGATGATGATT**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001205108
Insert Size: 513 bp
OTI Disclaimer:

Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

OTI Annotation:	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_001205108.1</u> , <u>NP_001192037.1</u>
RefSeq Size:	522 bp
RefSeq ORF:	513 bp
Locus ID:	668976
Cytogenetics:	X A1.1