

Product datasheet for **MC226058**

Gm3264 (NM_001242945) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Gm3264 (NM_001242945) Mouse Untagged Clone
Tag: Tag Free
Symbol: Gm3264
Synonyms: 100041306
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC226058 representing NM_001242945
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGACTAAGAAGAGATCAAAAAATAATGAACTAGAAGAACTGAAATTGGATATGAGGAAGATCAGCAATG
ACATGGAGGAAATGTGTGGAATCCTGAACGTTTACATGTATGAGGATTTGAACTACAGGATGAACACTGA
ATTCAACATCATTAAATCACAACATGAGAAGACAATGTTGGATATGAATAAAATGATCCAGTCCATAATT
GGTTCCATGCAGTACTCCAAGGAACTGATAGAAGATAACTATTCTACAGCATTAAAGGAGGACCCTCC
TCCGTGAGTGCCTCAACTCCACGAAAACGTAAGGATATTACTGAATGAGAACAGAAAGGCTGCTGGTGGG
GCAGGCTGGCCACAAGTGTCTGTGGGGAAGAAAAGAGTTTCTCTGAGGAGGCCAGCAAGAACATCTGTG
TCCAAGTGCCAAGGAACAGCAGTCCAGCAGAAAAGCAGAACATGGCACAGACCAGGACATGATCTCCCTC
AAAGAGAAGTGTGGAGGAAGAGCACTGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001242945
Insert Size: 519 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

OTI Annotation:	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001242945.1</u> , <u>NP_001229874.1</u>
RefSeq Size:	519 bp
RefSeq ORF:	519 bp
Locus ID:	100041306
Cytogenetics:	14