

Product datasheet for **MC226185**

Echdc2 (NM_001254754) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Echdc2 (NM_001254754) Mouse Untagged Clone
Tag: Tag Free
Symbol: Echdc2
Synonyms: 1300017C12Rik; 2610009M20Rik; D4Ertd765e
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC226185 representing NM_001254754
Red=Cloning site **Blue**=ORF **Orange**=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGAGTGATGTGGAGGTGGGACCTTTGTCCAGCGGCTCCGAGGCCTGATGAGTGAGATCGCAGCCTTCC
CTGTGCCACCATTGCAGCCATGGATGGGTTTGCCTTAGCGGGGGCTTGAACCTGCCCTGGCCTGTGA
CCTCCGAATAGCAGCTTCTTCGGCTGTCATGGGGCTAATTGAGACCACTCGAGGACTTCTCCAGGAGCA
GGAGGGACTCAGAGGCTGCCTCGCTGCCTGGGCGTGGCCCTAGCGAAGGAGCTCATCTTCACAGGCCGGC
GACTGAACGGGGCACAGGCCCGTGAGTTGGGCCTGGTAAATCATGCCGTGGCTCAGAATGAGGAGGGTAA
CGCTGCCTACCACGGGCCCTGGCACTGGCACAGGAGATCCTGCCCCAGGCTCCCATTGCTGTGCGCCTG
GGCAAAGTTGCCATCGACAGAGGAATGGAGGTGGACATTGCATCAGGGATGGCCATAGAACAGATGTGCT
ATGCCAGAACATCCCAACCCAGGACAGGCTGGAAGGAATGGCAGCCTTCAGGGAGAAACAGACTCCCAA
ATTTGTGGGCAAG**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: Sgfl-MluI
ACCN: NM_001254754
Insert Size: 576 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

| | |
|-------------------------------|---|
| OTI Annotation: | Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag. |
| Components: | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water). |
| Reconstitution Method: | <ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C. |
| RefSeq: | <u>NM_001254754.1</u> , <u>NP_001241683.1</u> |
| RefSeq Size: | 1255 bp |
| RefSeq ORF: | 576 bp |
| Locus ID: | 52430 |
| UniProt ID: | <u>Q3TLP5</u> |
| Cytogenetics: | 4 50.3 cM |