

Product datasheet for **MC226190**

Cyhr1 (NM_001276333) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Cyhr1 (NM_001276333) Mouse Untagged Clone
Tag: Tag Free
Symbol: Cyhr1
Synonyms: 1110031M01Rik; AU042374; Chrp
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC226190 representing NM_001276333
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGTTTCCAACCAGGACTGAATGGAGCACTGCTGTCCACCTGGTATTGGCAGGGGTGTCCTGC
ATGCAGCGGTGAGCTCTGTACAGTCCAGTCGAGGAGCAGCTGCCGGCTTCTGCTCCAGACCTTTGCTGC
CATCATTATGTTGGCCCCAGGGCCGAGCACACATGAAGACTGTCTCGCTGGAGCTTGGGTTGCCACAGTC
ATCGGTCTGCCTTCTGGCCTTCGATTTTCACTGGGTGAATGGGACCGCTTCTGCCAACCTGCCTC
TGGGAGGAGGCATGGTGCTGGCAGTGGCTGGTGATCACCTGGGCCCTGAAGGGTGCTCTGTGCTGGGCA
GGCAGTACTGCTGGTGGTTGCAGTGACCATCCTTATTGTGGCTGTTTTACAGCTAACACTTACGGGATG
TGGGGGGGTACCTGCTGGGTGTGGCAGGCCTTGTAGCAGGCTGGAAGAAGACAGGCTGCTACTGCTGC
CAAAGGAGGATGTCTGTGCTGGGCCCTGGCTGCAGGTAGTTGGGCTATTGCAGGGCCCTGCACACACA
GCGCCTGCAGTGGGAG**TGA**

ACGCGTACGCGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001276333
Insert Size: 579 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

OTI Annotation:	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001276333.1</u> , <u>NP_001263262.1</u>
RefSeq Size:	1267 bp
RefSeq ORF:	579 bp
Locus ID:	54151
UniProt ID:	<u>Q9QXA1</u>
Cytogenetics:	15 D3