

## Product datasheet for **MC226191**

### Cyhr1 (NM\_001276334) Mouse Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Cyhr1 (NM_001276334) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Cyhr1
Synonyms:	1110031M01Rik; AU042374; Chrp
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
Fully Sequenced ORF:	>MC226191 representing NM_001276334 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGTTTCCAACCAGGACTGAATGGAGCACTGCTGTCCACCTGGTATTGGCAGGGGTGTCCTGC  
ATGCAGCGGTGAGCTCTGTACAGTCCAGTCGAGGAGCAGCTGCCGGCTTCTGCTCCAGACCTTTGCTGC  
CATCATTATGTTGGCCCCAGGGCCGAGCACACATGAAGACTGTCTCGCTGGAGCTTGGGTTGCCACAGTC  
ATCGGTCTGCCTTCTGGCCTTCGATTTTCACTGGGTGAATGGGACCGCTTCTGCCAACCTGCCTC  
TGGGAGGAGGCATGGTGCTGGCAGTGGTGGTGATCACCTGGGCCTGAAGGGTGCTGTGCTGCTGGGCA  
GGCAGTACTGCTGGTGGTTGCAGTGACCATCCTTATTGTGGCTGTTTTACAGCTAACACTTACGGGATG  
TGGGGGGGTACCTGTGGTGTGGCAGGCCTTGTAGCAGGCTGGAAGAAGACAGGCTGCTACTGCTGC  
CAAAGGAGGATGTCTGTCTGGTGGCCCTGGCTGCAGGTAGTTGGGCTATTGCAGGGCCCTGCACACACA  
GCGCCTGCAGTGGGAG**TGA**

**ACGCGT**ACGCGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	Sgfl-Mlul
ACCN:	NM_001276334
Insert Size:	579 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

<b>OTI Annotation:</b>	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001276334.1</a></u> , <u><a href="#">NP_001263263.1</a></u>
<b>RefSeq Size:</b>	1051 bp
<b>RefSeq ORF:</b>	579 bp
<b>Locus ID:</b>	54151
<b>UniProt ID:</b>	<u><a href="#">Q9QXA1</a></u>
<b>Cytogenetics:</b>	15 D3