

Product datasheet for **MC226231**

Tk1 (NM_001271729) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Tk1 (NM_001271729) Mouse Untagged Clone
Tag: Tag Free
Symbol: Tk1
Synonyms: D530002A18Rik; Tk-1; Tk1a; Tk1b
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC226231 representing NM_001271729
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAGAAGAGTCCGGCGCTTCCAGATCGCCAGTACAAGTGCCTGGTCATCAAGTATGCCAAGACACGC
GCTATAGCAACAGCTTCCACACATGATCGGAACACCATGGACGCATTGCCAGCCTGCATGCTCCGCGA
TGTGACCCAGGAGTCCTGGGTGTGGCCGTCATTGGCATCGATGAGGGGAGTCTTTTCTGACATTGTG
GATTTCTGTGAAATGATGGCCAACGAGGGCAAGACAGTAATTGTGGCAGCGCTGGATGGGACCTCCAGA
GGAAGGCTTTCGGCAGCATCTTGAACCTGGTGCCTGGCGGAGAGTGTGGTGAAGCTCACCGCTGTGTG
CATGGAGTGCTCCGAGAAGCTGCCTACACGAAGAGGCTGGGCTGGAGAAAGAGGTGGAGGTGATTGGC
GGAGCCGACAAGTATCACTCCGTGTGCCGCTGTGCTACTTTAAGAAGTCTTCAGCCAGACTGCTGGCT
CAGACAACAAGAACTGTCTGGTGTGGGCGAGCCGGGAGAGGCCCTGGTTGTGAGGAAGCTTTGCCTC
TCAGCAAGTCCTACAATACTGCAACT**TGA**

ACGCGTACGCGCGCCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001271729
Insert Size: 594 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

OTI Annotation:	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001271729.1</u> , <u>NP_001258658.1</u>
RefSeq Size:	1357 bp
RefSeq ORF:	594 bp
Locus ID:	21877
UniProt ID:	<u>P04184</u>
Cytogenetics:	11 82.96 cM