

Product datasheet for **MC226257**

Arhgdib (NM_001301303) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Arhgdib (NM_001301303) Mouse Untagged Clone
Tag: Tag Free
Symbol: Arhgdib
Synonyms: D4; Gdid; Gdid4; Ly-G; Ly-GDI
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC226257 representing NM_001301303
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGACGGAGAAGGATGCACAGCCACAGCTGGAAGAGGCGGACGACACCTGGACAGCAAGCTCAATTATA
AGCCACCCCTCAGAAGTCCTTGAAGGAGCTGCAGGAGATGGACAAGGATGACGAGAGTCTAACCAAGTA
CAAGAAAACACTGCTGGGAGATGTCCTGTGGTAGCAGACCCAACAGTCCCAACGTGACTGTTACCCGG
CTTAGCCTTGATGTGACAGTGACACAGGCCATCACCATGGACCTTACTGGCGATCTCGAGGCCCTCA
AAAAGGATACATTTGTGCTAAAGGAAGGCATTGAATACAGGGTAAAAATTAACCTCAAAGTGAATAAGGA
TATTGTGTCTGGCCTGAAGTATGTTCAACACACATACCGGACTGGCATGAGAGTGGATAAAGCCACATTC
ATGTTGGCAGCTATGGGCCCGACAGAGGAGTACGAATTCCTCACTCCAGTAGAGGAAGCTCCCAAGG
GCATGCTGGCCCGAGGCACTTACCACAACAAGTCCTTCTCACGGATGACGACAAACAGGACCACCTCAC
CTGGGAATGGAACCTGGCCATTAAGAAGGATTGGACAGAA**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001301303
Insert Size: 603 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

OTI Annotation:	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001301303.1 , NP_001288232.1
RefSeq Size:	1170 bp
RefSeq ORF:	603 bp
Locus ID:	11857
UniProt ID:	Q61599
Cytogenetics:	6 G1
Gene Summary:	<p>The protein encoded by this gene is a member of the Rho guanine nucleotide dissociation inhibitor (GDI) family. This gene is expressed at high levels in hematopoietic cells. This protein is cytosolic, and dissociation of Rho from this protein is required for membrane association and activation of Rho by Guanine Nucleotide Exchange Factors (GEFs). C-terminal truncations of this gene product have been reported to promote metastasis. Multiple transcript variants and protein isoforms exist. [provided by RefSeq, Aug 2014]</p> <p>Transcript Variant: This variant (3) contains an alternate exon in the 5' UTR compared to variant 1. Variants 1, 2, 3 and 4 encode the longer protein, isoform 1. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>