

Product datasheet for MC226305

Nudt7 (NM_001290180) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Nudt7 (NM_001290180) Mouse Untagged Clone
Tag: Tag Free
Symbol: Nudt7
Synonyms: 1300007B24Rik; 2210404C19Rik
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC226305 representing NM_001290180
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTGAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGC**

ATGCGCGTCAGAACTATTTGCTTCCGAAGTCTCCTGCTGCTAGTAAGAGGAAGTACTGCAGGGCTGATG
 TCGCGACCTTGTGACTCCCGAGCCTGTCAGCTGAAAAGGGAACCTGGAGAAGTCTGCTTCCCTGGAGG
 AAAGAGGGACCCGGTGGACACAGATGACACAGCCACTGCTCTCCGTGAAGCCAGGAGGAGTGGGGCTG
 CATCCCCACCAAGTGGAGGTGGTCTCTACCTGGTGCCATACGTATTTGATAATGATGCACTGGTAACCC
 CCGTAGTGGGTTTTCTAGACCACAACTCCAGGCCCAACCTAATGCTGATGAAGTGAAGGAAGTCTTCTT
 TGTGCCTTTGACTATTTCTCCATCCCAAGTCTACTACCAGAAGCAAATCACACAGTCTGGCCGTGAT
 TTCATCATGCATTGCTTCGAGTACAAAGACCCTGAGACTGGTGTGAACCTACCTAATCCAGGGAATGACCT
 CAAACTGGCTGTGTGGTGGCCTTAATTATTTGGAACAAAGTCTGCCTTCAAGATTGATTTTGATCT
 CCATGACCTGATACCGTCTTGTAAGGACCTTCCTTTGAGATATTCTTTAAGCAAGTTG**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001290180
Insert Size: 624 bp
OTI Disclaimer:

Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).


[View online »](#)

OTI Annotation:	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001290180.1</u> , <u>NP_001277109.1</u>
RefSeq Size:	945 bp
RefSeq ORF:	624 bp
Locus ID:	67528
UniProt ID:	<u>Q99P30</u>
Cytogenetics:	8 E1
Gene Summary:	<p>Coenzyme A diphosphatase which mediates the cleavage of CoA, CoA esters and oxidized CoA with similar efficiencies, yielding 3',5'-ADP and the corresponding 4'-phosphopantetheine derivative as products. CoA into 3',5'-ADP and 4'-phosphopantetheine. Has no activity toward NDP-sugars, CDP-alcohols, (deoxy)nucleoside 5'-triphosphates, nucleoside 5'-di or monophosphates, diadenosine polyphosphates, NAD, NADH, NADP, NADPH or thymidine-5'-monophospho-p-nitrophenyl ester. May be required to eliminate oxidized CoA from peroxisomes, or regulate CoA and acyl-CoA levels in this organelle in response to metabolic demand. Does not play a role in U8 snoRNA decapping activity. Binds U8 snoRNA. [UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (3) lacks an alternate exon in the 5' coding region, resulting in translation initiation at an alternate start codon, compared to variant 1. The encoded isoform (3) has a distinct N-terminus and is shorter than isoform 1.</p>