

Product datasheet for **MC226323**

Ahsg (NM_001276449) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Ahsg (NM_001276449) Mouse Untagged Clone
Tag: Tag Free
Symbol: Ahsg
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC226323 representing NM_001276449
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGAAGTCCCTGGTCTTGCTCCTTTGTTTTGCTCAGCTCTGGGGCTGCCAATCCGCTCCACAAGGTACAG
GACTGGGTTTTAGAGAATTGGCTTGTGATGATCCAGAAGCAGAGCAAGTAGCTTTGTTGGCCGTGGACTA
CCTCAATAATCATCTTTCAGGGATTCAAACAGGTCTTGAATCAGATCGACAAAGTCAAGGTGTGGTCT
CGGCGGCCCTTCGGAGTGGTGTATGAGATGGAAGTTGACACACTGGAGACCACTTGCCATGCTTTGGACC
CCACCCGCTGGCAAATCTCATGCATAATCTTGGTGGGGAAGAAGTTTCAGTGGCCTGCAAGTTATTCCA
AACACAGCCCCAGCCAGCCAATGCCAACGCAGTAGGTCCCGTACCCACAGCGAATGCAGCCCTACCAGCT
GACCCACCTGCATCTGTGGTGGTGGGACCTGTGGTGGTCCACGAGGACTTTCAGACCACCGAACTTACC
ACGACCTACGCCACGCCTTCTCTCTGTGGCCTCGGTGGAGTCCGGCTCGGGAGAACTTTCATTCTCC
TAAGGTGGCCAGCCTGGTGTCTGTGGTCCAGTGTCCCCATGTGCCAGGGAGGATCAGACACTTCAA
ATCTAG

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001276449
Insert Size: 636 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

OTI Annotation:	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001276449.1</u> , <u>NP_001263378.1</u>
RefSeq Size:	1122 bp
RefSeq ORF:	636 bp
Locus ID:	11625
Cytogenetics:	16 13.74 cM
Gene Summary:	Probably involved in differentiation.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (2) has multiple differences in the coding region but maintains the reading frame, compared to variant 1. This variant encodes isoform 2, which is shorter than isoform 1.