

Product datasheet for **MC226364**

Sigmar1 (NM_001286538) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Sigmar1 (NM_001286538) Mouse Untagged Clone
Tag: Tag Free
Symbol: Sigmar1
Synonyms: mSig; O; Oprs1; Si; Sig1R; sigma1R
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC226364 representing NM_001286538
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGCCGTGGGCCGCGGGACGGCGGTGGGCATGGATCACCTGATTCTGACTATTATCGCAGTGCTGATCC
AGGCCGCTGGTTGTGGCTGGGCACTCAAACCTTCGTCTTCTCTAGAGAAGAAATAGCGCAGCTTGCTCG
ACAGTATCGGGGCTGGACCATGAGCTTGCTTCTCTCGGCTGATCGTGGAGCTGCGGAGGCTGCACCCA
GGCCACGTGCTGCCGGATGAGGAGCTGCAGTGGGTATTTGTGAACGCGGGCGGCTGGATGGCGCCATGT
GTATTCTGCACGCCTCGCTGTCTGAGTACGTGCTGCTCTTCGGCACCGCCCTGGGCTCCCATGGCCATTC
GGGACGATACTGGGCTGAGATTTCTGACACCATCATCTCTGGCACCTTCCACCAATGAAAGAGGGCACC
ACGAAAAGTGAGGTCTTCTACCCAGCACCCAGGACTACCTCACACTCTTCTATACCCTTCGGGCCTATGC
CCGGGGCCTCCGGCTTGAGCTTACCACCTACCTCTTTGGCCAAGACTCCTGACCAGCCAGGCTGAAGGA
AGATCTGTGGATGGACAGGAGCGCGGACGAGCACATCCATCTCTCTACTCGCTGGAACCTCATATTCAGA
CAGCAGCGTATACCATGCAGATAT**TAA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001286538
Insert Size: 657 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001286538.1</u> , <u>NP_001273467.1</u>
RefSeq Size:	1510 bp
RefSeq ORF:	657 bp
Locus ID:	18391
UniProt ID:	<u>O55242</u>
Cytogenetics:	4 A5
Gene Summary:	<p>This gene encodes a transmembrane protein located in the endoplasmic reticulum. The encoded protein is a receptor that binds several endogenous ligands, including N,N-dimethyltryptamine, progesterone and pregnenolone and a variety of of non-opiate compounds. The encoded protein plays a role in regulating the activity of ion channels, acting as a chaperone and protecting cells from oxidative stress. In humans, this receptor has been associated with Alzheimer's and Parkinson's diseases, stroke and numerous disease conditions such as depression, pain and addiction. Alternative splicing results in multiple transcript variants encoding different isoforms.[provided by RefSeq, Nov 2013]</p> <p>Transcript Variant: This variant (2) differs in the 3' UTR and uses an alternate splice site in the 3' terminal exon which results in a frameshift compared to variant 1. The encoded protein (isoform 2) is shorter and has a distinct C-terminus compared to isoform 1.</p>