

Product datasheet for **MC226365**

Sdf4 (NM_001302469) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Sdf4 (NM_001302469) Mouse Untagged Clone
Tag: Tag Free
Symbol: Sdf4
Synonyms: Cab4; Cab45
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC226365 representing NM_001302469
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGTCTGGCTGGTGGCAATGACGCCAGGCAGAGTTCCTCTGTGGTCTAGCTGCTCATGGCCTCTGGT
 TCTTGGGCCTTGTCTTCTGATGGATGCAACCGCTAGACCTGCCAACCACTCATCTACTCGGAAAGAGC
 GGCTAACAGGGAGGAAAATGAGATCATGCCCCAGACCACCTGAATGGGGTGAAGCTGGAGATGGATGGA
 CACCTCAACAAGGACTTCCATCAGGAGTTTTCTGGGAAAGGACATGGATGGGTTTGATGAAGACTCAG
 AGCCACGGAGAAGCCGGAGGAAGCTGATGGTCATCTTTCCAAGGTAGACGTGAACACTGACCGGAGGAT
 CAGCGCTAAGGAGATGCAGCACTGGATTATGGAGAAAACAGCAGAGCACTTCCAGGAGGCCGTCAAGGAA
 AACAACTGCACTTCAAGGCTGTGGACCCTGACGGTGACGGCCATGTTTCTGGGATGAATATAAAGTGA
 AGTTTTTGGCAAGCAAAGGCCACAATGAAAGGGAGATTGCTGAAGCCATCAAGAACCATGAGGAGCTCAA
 AGTGGATGAGGAGATTGCGCTTCATTTCTGCTGCTTACTGGAGGGCCCTGTCTGTCTCTGGTGAG
 ACCAGTATATGCCTGTCCAGTTG**TAG**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001302469
Insert Size: 657 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001302469.1 , NP_001289398.1
RefSeq Size:	990 bp
RefSeq ORF:	657 bp
Locus ID:	20318
UniProt ID:	Q61112
Cytogenetics:	4 E2
Gene Summary:	<p>This gene encodes a member of the CREC family. The encoded protein contains multiple calcium-binding EF-hand motifs. This protein localizes to the Golgi lumen and may be involved in regulating calcium dependent cellular activities. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Oct 2014]</p> <p>Transcript Variant: This variant (4) uses an upstream alternate 3' terminal exon compared to variant 1. It encodes isoform 2, which has a shorter and distinct C-terminus, compared to isoform 1.</p>