

## Product datasheet for **MC226368**

### Gimap4 (NM\_001243199) Mouse Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Gimap4 (NM\_001243199) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Gimap4  
**Synonyms:** AU019574; E430007K16Rik; I; lan-1; lan1; lmap4  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >MC226368 representing NM\_001243199  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGAAGTCCAGTCCGGTGGTGC GGTTTCATCCAGAAAGTTCAAGGAGCAGCCATGAGCTTGAAACC  
AAGATCAAGGAATCCCCAACTGAGAATTGTCTTACTTGGAAAACTGGAGCAGGAAAGAGTTCAACAGG  
GAACAGTATCCTTGGGAAAAAGTGTAAATTCTGGCATTGTGCAAAATCCATCACCAAGTCTGTGAA  
AAAAGGTGAGCACCTGGGATGGGAAAGAGCTTGTGTCGTGGATACACCTGGTATTTTGCACACTGAGG  
TACCAGATGCTGACACACAAAGGAGATCACTCGCTATGTTGCCCTGACCTCTCCAGGGCCTCATGCTCT  
GCTCCTGGTAGTCCACTGGGGCGTTACTGTGGAAGAACAAGGCTACACAGAAAATCTGGACATG  
TTTGGAAACAGGCTAGAAGATTCATGATTCTCTTGTCTCACCAGGAAGGATGACTTAGAAGACACTGATA  
TCCATGAGTACTTAGAGAAGGCTCCTAAATCTTTCAAGAGGTGATGCATGAGTCCAGAATCGCTACTG  
TTTGTTCACAACAGAGCCTCAGGCGTGAAAAGGAAGAGCAGAAGATGCAATTGTTGACCTTGGTCCAG  
AGCATGTTTCTGTCTCCAGGATGAAG**TAA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_001243199  
**Insert Size:** 660 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<u><a href="#">NM_001243199.1</a></u> , <u><a href="#">NP_001230128.1</a></u>
<b>RefSeq Size:</b>	1365 bp
<b>RefSeq ORF:</b>	660 bp
<b>Locus ID:</b>	107526
<b>UniProt ID:</b>	<u><a href="#">Q99JY3</a></u>
<b>Cytogenetics:</b>	6 B2.3
<b>Gene Summary:</b>	<p>This gene encodes a protein belonging to the GTP-binding superfamily and to the immuno-associated nucleotide (IAN) subfamily of nucleotide-binding proteins. This gene exists within a cluster of other related genes located on mouse chromosome 6. This family member encodes a lymphoid signaling protein that functions to accelerate programmed T-cell death, which appears to correlate with the phosphorylation status of the protein. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Aug 2011]</p> <p>Transcript Variant: This variant (3) lacks an alternate segment that results in a frameshift in the 3' coding region, compared to variant 1. The encoded isoform (c) has a distinct C-terminus and is shorter than isoform a. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>