

Product datasheet for **MC226377**

TPD52I2 (NM_001291197) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Tpd52I2 (NM_001291197) Mouse Untagged Clone
Tag: Tag Free
Symbol: Tpd52I2
Synonyms: 2810411G23Rik; AU016537; C86069; D54
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC226377 representing NM_001291197
Red=Cloning site **Blue**=ORF **Orange**=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGGACTCTGCTAGCCAAGATATCAACCTGAATTCCTCAACAAAGGTGTGCTGTCTGACTTTATGACTG
ACGTCCCTGTTGACCCAGGTGTGGTCCACCGACTCCAGTTGTAGAAGGCCTGACAGAGGGGGAGGAAGA
AGAGCTTCGGGCTGAGCTTGCTAAGGTGGAAGAAGAAATTGCTACTCTGCGCCAGGTGCTGGCAGCCAAA
GAGAGGCACTGTGGAGAGCTGAAAAGGAGGCTGGGCCTCTCCACATTAGGGGAGCTGAAGCAGAACCTGT
CTAGGAGCTGGCATGATGTGCAGGTCTCTACTGCCTATGTGAAAACGCTCTGAGAACTTGAGAGTGGAA
TGAGAAAGTGACGCAGTCTGACCTCTACAAGAAGACTCAAGAACTCTTTCACAGGCTGGACAGAAAACA
TCAGCTGCCCTGTCCACCATGGGCTCTGCCATCAGCAGGAAGCTTGAGACATGAGCAGCAGCTACTCCA
TCCGCCACTCGATAAGTATGCCTGTGATGAGGAAGCTCAGCCACCTTCAAGTCATTTGAGACCGAGTGGG
GACCATAAAGTCTAAGTTGTGGGTGGCAGAGAGAATGGCAGCGATAACCTCCCTCCCTCTCTGGAAGT
GGTGACCAGACATTGCCGGATCATGCGCCTTTC**TAA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001291197
Insert Size: 666 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001291197.1</u> , <u>NP_001278126.1</u>
RefSeq Size:	3625 bp
RefSeq ORF:	666 bp
Locus ID:	66314
UniProt ID:	<u>Q9CYZ2</u>
Cytogenetics:	2 H4