

Product datasheet for **MC226407**

Spc25 (NM_001199123) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Spc25 (NM_001199123) Mouse Untagged Clone
Tag: Tag Free
Symbol: Spc25
Synonyms: 2600017H08Rik; 2610205L13Rik; Spbc; Spbc25
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC226407 representing NM_001199123
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGGGTGAAGACGAATTGGCACTCTTAATCAAAGCATCAATGAATTTGGGGATAAGTTCAGAAACAGAC
TCGATGATAATCATAGTCAAGTGTGGGACTTAGAGACGCCTTCAAGGACTCCATGAAAGCGTTTTTCAGA
AAAGATGTCTTTGAAATTAAGGAAGAAGAGAGGATGACTGAGATGATTCTGGAGTATAAAAACAGCTC
TGTAAGCAGAATAAGCTCATTCAAGAAAAGAAAGAGAATGTGTTGAAGATGATTGCTGAAGTAAAAGGCA
AGGAGCAAGAGTCGGAAGAGCTGACTGCTAAAATCCAGGAGCTCAAGGAAGAGTACGCTAGGAAGAGGGA
AACCATTTCCACTGCTAACAAAGCTAATGAAGAGAGATTGAAAGGACTGCAGAAATCAGCGGATCTGTAT
AGAGATTACCTTGGACTAGAAAATTAGAAAAGATTCACGGAATAAAATGCAGTTTTATTTACTAGTATTG
ACCCTAAGAAATCCTGAGAGCCCATATATGTTTTCCATGAGCATAAATGAAGCTAAGGAATATGAAGTGTA
CGACAGTTCGCCCTCATCTTGAGTGCCTAGCAGAATTTCCAGGAGAAAGTCAAGGAAGACCAACAATTTTCA
GCTTTTCTTGCCAATATTCGGAAGGCTTTTATAGCTAAGGTTTATAAT**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001199123
Insert Size: 681 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001199123.2</u> , <u>NP_001186052.1</u>
RefSeq Size:	1622 bp
RefSeq ORF:	681 bp
Locus ID:	66442
UniProt ID:	<u>Q3UA16</u>
Cytogenetics:	2 C2
Gene Summary:	<p>This gene encodes a component of the kinetochore-associated NDC80 protein complex, which is required for the mitotic spindle checkpoint and for microtubule-kinetochore attachment. During meiosis in mouse, the protein localizes to the germinal vesicle and then is associated with the chromosomes following germinal vesicle breakdown. Knockdown of this gene in oocytes results in precocious polar body extrusion, chromosome misalignment and aberrant spindle formation. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Apr 2015]</p> <p>Transcript Variant: This variant (1) represents the longest transcript and encodes the longer isoform (1). Variants 1, 2, and 3 encode the same isoform.</p>