

## Product datasheet for **MC226436**

### Il15ra (NM\_001271499) Mouse Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Il15ra (NM\_001271499) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Il15ra  
**Synonyms:** AA690181; IL-15RA  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >MC226436 representing NM\_001271499  
**Red**=Cloning site **Blue**=ORF **Orange**=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**GCGATCGCC**

ATGGCCTCGCCGAGCTCCGGGGCTATGGAGTCCAGGCCATTCTGTGTTGCTGCTGCTGCTGTTGCTAC  
TGTTGCTCCCGCTGAGGGTGACGCCGGCACCACGTGTCCACCTCCCGTATCTATTGAGCATGCTGACAT  
CCGGGTCAAGAATTACAGTGTGAACTCCAGGGAGAGGTATGTCTGTAACCTGGCTTTAAGCGGAAAGCT  
GGAACATCCACCCTGATTGAGTGTGTGATCAACAAGAACAACAAATGTTGCCCACTGGACAACCTCCAGCC  
TCAAGTGCATCAGAGACCCCTCCCTAGCTCACTACAGTCCAGTGCCAACAGTAGTGACACCAAAGGTGAC  
CTCACAGCCAGAGAGCCCTCCCTCTGCAAAAGGAACTACAGGGACAGGCAGTCACAAGTCTCCCGA  
GCCCATCTCTTGAGCAACAATGACCTTGAGCCTACAGCCTCCACCTCCCTCAGGATAACAGAGATTT  
CTCCCCACAGTTCCAAAATGACGAAAGTGGCCATCTCTACATCGGTCCTTTGGTTGGTGCAGGGGTTGT  
GATGGCTTTCTGGCCTGGTACATCAAATCAAGGCAGCCTTCTCAGCCGTGCCGTGTTGAGGTGAAACC  
ATGGAACAGTACCAATGACTGTGAGGGCCAGCAGCAAGGAGGATGAAGACACAGGAGCC**TAA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_001271499  
**Insert Size:** 693 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001271499.1</a> , <a href="#">NP_001258428.1</a>
<b>RefSeq Size:</b>	1565 bp
<b>RefSeq ORF:</b>	693 bp
<b>Locus ID:</b>	16169
<b>UniProt ID:</b>	<a href="#">Q60819</a>
<b>Cytogenetics:</b>	2 8.97 cM
<b>Gene Summary:</b>	<p>High-affinity receptor for interleukin-15 (PubMed:17947230). Can signal both in cis and trans where IL15R from one subset of cells presents IL15 to neighboring IL2RG-expressing cells (PubMed:17947230). In neutrophils, binds and activates kinase SYK in response to IL15 stimulation (By similarity). In neutrophils, required for IL15-induced phagocytosis in a SYK-dependent manner (By similarity).[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (5) uses an alternate in-frame splice site in the coding region compared to variant 1. It encodes isoform 4 (also known as isoform 1A) which is shorter than isoform 1. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>