

Product datasheet for **MC226448**

Ube2j2 (NM_001284314) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Ube2j2 (NM_001284314) Mouse Untagged Clone
Tag: Tag Free
Symbol: Ube2j2
Synonyms: 1200007B18Rik; 2400008G19Rik; 5730472G04Rik; AL022923; NCUBE-2; Ubc6; Ubc6p
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC226448 representing NM_001284314
Red=Cloning site **Blue**=ORF **Orange**=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGAGCAATAACAGCAATAAGAGAGCTCCAACAACAGCAACACAAAGGTTGAAGCAGGACTACCTTCGAA
TTAAGAAAAGCCAGTGCCTTACATCTGTGCAGAGCCCTCCCTTCAAATATTCTTGAATGCATTATGT
TGTCAGGAGCCCTGAGATGACTCCTTATGAAGGTGGCTATTATCATGGAAAATAATTTTCCAGAGAA
TTTCCATTTAAGCCTCCTAGTATTTACATGATAACCCCAATGGAAGATTTAAGTGCAACACAAGGCTGT
GTCTCTCCATCACAGATTTCCACCCAGATACATGGAATCCAGCGTGGTCCGTCTCCACCATTCTAACAGG
GCTCCTGAGCTTTATGGTGGAGAAAGGCCCCACCTGGGGAGCATAGAGACGTCGGACTTCACGGAAATT
AAACAGAAACAGAAAGCACAAGATGAACAAAGTAAACAGGCCCCAGAACTTTCCTTTGCCAGACGTGGTTC
CTGATGGGGAGCTGCACCGTGGCCAGCATGGGATCCAGCTCCTCAATGGGCACGCACCCGCTGTGGCC
CAACCTCGCTGGGCTCCACAGGCCAACAGGCATCACGGACTCCTGGCGGCGGCTTGGCGAACTTGTTT
GTTATAGTTGGGTTTGCAGCCTTCGCCTACACGGTCAAGTACGTGCTGAGAAGCATAGCACAGGAG**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001284314
Insert Size: 699 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001284314.1</u> , <u>NP_001271243.1</u>
RefSeq Size:	3139 bp
RefSeq ORF:	699 bp
Locus ID:	140499
UniProt ID:	<u>Q6P073</u>
Cytogenetics:	4 87.66 cM
Gene Summary:	<p>Catalyzes the covalent attachment of ubiquitin to other proteins. Seems to function in the selective degradation of misfolded membrane proteins from the endoplasmic reticulum (ERAD).[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (6) differs in the 5' UTR, 5' and 3' coding regions and initiates translation at a downstream start codon, compared to variant 1. It encodes isoform c, which is shorter than isoform a. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>