

Product datasheet for **MC226449**

Serpina3a (NM_028740) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Serpina3a (NM_028740) Mouse Untagged Clone
Tag: Tag Free
Symbol: Serpina3a
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC226449 representing NM_028740
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAGGACAGTCTGGCTATTTAGATGTTTCCTTTTTGCTGGGTCCTCAATATCAGACAGGAGCTACTTG
AAGGCTTTTGAACGTGACTTTTGACCCTGAAGACACGTTCTGGGAAATTTACATTGGATAGGAAGAG
GACCGTGAATGTGCCATGATGAAAACCGAGGAGCTAACCAAACTACTCCGGGACGAGGAGATGCAA
AGTACTGTGATGGAGCTGAATTACATAGGCAATGCCAGCTTCTGTTTATCCTCCCTGACCAGGGCAGGA
TACAGCACGTGGAAGATAGCTTGAACCCAGTCTCTGAGGAAATGGAGGAAATCTTTGAGGCCAGGAT
GCTAGATGAACTCTCCCTGCCAAGTTCTCCTTATCACAAGACTACAACCTAAATGACATCCTTCCAGAG
CTGGGTATCAAGGAAGTCTTCTCCACACAGGCTGACCTGTCTGGGATCACTGGAGCCAAAAATAAAGAG
TCTCTCAGATGATCCACCAGGCTGCGCTGGATGTTACCGAGACACACAGAAGCAGATGTCATCACAAT
AGCCCGATATAACTTCCAATCTGCAAAAATTAAGGCTAAAATTGTTAAGGTGGACAGGGAATTTGTAC
CTTATTTTAGACCAATGTTAAGAGTATTAGCGTTATGGGTAAGGTTATCAACCCTTAAACAACTAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_028740
Insert Size: 699 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

OTI Annotation:	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_028740.2</u> , <u>NP_083016.1</u>
RefSeq Size:	1204 bp
RefSeq ORF:	699 bp
Locus ID:	74069
UniProt ID:	<u>Q6P4P1</u>
Cytogenetics:	12 E