

Product datasheet for **MC226451**

Stard6 (NM_001289648) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Stard6 (NM_001289648) Mouse Untagged Clone
Tag: Tag Free
Symbol: Stard6
Synonyms: 1700011K09Rik; 4833424I06Rik; 4933429L05Rik; AI644424
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC226451 representing NM_001289648
Red=Cloning site **Blue**=ORF **Orange**=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGACTATAAGGCAATTGCCCAACAACTGCTGAACAAGTTTTAGCTTACAATCAAGACCTATCTGGCT
GGAAACTCATTAAATCTCAAAAAAGGTAAGTGTTCAGCAAGACCTCTAGAATATCCACGGAAACCT
ATACCGCGTTGAAGGAATAATCCAGAATCAGCAGCTCATCTCTGACTTCCTCTCAAACATGATCAC
AGAGTCTCATGGGACAAATCATTGAAAGGTTCAATGTGATACACAAGATTGATTCGGACACACTCATAT
GCCACACCATTACACAAAGTTTTGCCATGGGCTCAATTTCCCCTAGAGACTTTATCGACTTAGTACACAT
CAAGCACTATGAAAGGAACGTGGACATTATCAGTACTAAAAGTGTGGACTTCCCAGGATATGCTCCAAC
TCGACTTATATCCGTGGCTTTAATCATCCCTCTGGCTATGTGTGTTACCTCTGAAGGAAAACCCAGCCT
ATTCCAAGCTAGTGATTTTTGTCCAGACAGAGATGAAAGGAAAATTGCCTGCATCGGTAATTGAAAAGTC
TATGCCTTCCAACCTAGTAAGCTTCTCCTCAATGTGAAAGATGGCGTAAAGACGTACAGAATCCACCC
ATACGTGCCCGGCACAGTTCGACTCCTCTGTGCACAAAAAGAAAGGACATTCTGCCATTAACCGT
AG

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001289648
Insert Size: 702 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001289648.1</u> , <u>NP_001276577.1</u>
RefSeq Size:	1427 bp
RefSeq ORF:	702 bp
Locus ID:	170461
UniProt ID:	<u>P59096</u>
Cytogenetics:	18 44.44 cM
Gene Summary:	May be involved in the intracellular transport of sterols or other lipids. May bind cholesterol or other sterols (By similarity).[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (2) differs in the 5' UTR, compared to variant 1. Variants 1, 2, 3, and 4 encode the same protein.