

Product datasheet for **MC226461**

Lin37 (NM_001290569) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Lin37 (NM_001290569) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Lin37
Synonyms:	1810054G18Rik
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>MC226461 representing NM_001290569 Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**

ATGGAGATGGCCAAAGCCCGGAACCCAGCTGGATGCCGTTCTTCAGTGTCTGCTGGAGAAGAGTCACATGG
ACAGGGAACGCCTGGATGAAGAAGCTGGGAAAACCCCTTGGATACGCACAACAAGGATTGTTCCATTGC
TGCCACAGGCAAAAGGCCATCTGCCGCTTCCCCACGCGGAGGAAGAAGAGGAGGGAGATGGACGAT
GGGCTCGCTGAGGGGGTCTCAGAGATCCAATACGTATGTGATCAAGCTATTTGACCGAGTGTGGACT
TGGCCAGTTCAGTGAGAACACACCACTGTACCCCATCTGCCGCGCTGGATGCGCAACAGCCCCACAGT
GCGAGAGCGCGAGCGCTACCCGGCTACCGCTGCCTCCTCTGCCGAGGATGGAGAGGGTTCAGAGGTC
ATCAACAGCAAAAACCGTGATGTGTACAAGCTGCCTCCACCCACAGCTCCTGGACCACTTGAGATGCCT
GTAGATCCAGAATTCATCCCCACTGCAGCCTGAGACCGAGGGTACCCCGATGATGAACCTTCCGAGCC
TGAACCTTCTCCTTCCACACTTATCTATCGCAACATGCAACGTTGGAAACGCATACGCCAGAGGTGGAAG
GAGGCATCTCATCGGAACCAGCTTCGCTACTCAGAAAGTATGAAGATCCTACGGGAGATGTATGACCGGC
AGTGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	SgfI-MluI
ACCN:	NM_001290569
Insert Size:	705 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_001290569.1, NP_001277498.1</u>
RefSeq Size:	1124 bp
RefSeq ORF:	705 bp
Locus ID:	75660
UniProt ID:	<u>Q9D8N6</u>
Cytogenetics:	7 B1