

Product datasheet for MC226579

Kctd17 (NM_001289671) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Kctd17 (NM_001289671) Mouse Untagged Clone
Tag: Tag Free
Symbol: Kctd17
Synonyms: 2900008M13Rik; AA414907; AW742389; N28155
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC226579 representing NM_001289671
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**GCGATCGCC**

ATGCAGACAACCGCGCCGGCGATGAGGATGGAGGCCGGGAGGCAGCGCCACCGGTGGGGCGGGCGGCC
 GCCCGGGGGCGGCTGGGGCAAGTGGGTGCGGCTGAACGTGGGGGCACGGTGTTCCTGACCACCCGCCA
 GACGCTGTGTCGGGAGCAGAAGTCTTCCTCAGCCGCTGTGTCAGGGCGAAGAGCTTCAATCGGACAGG
 GATGAGACCGGGGCTACCTCATTGACCGTGACCCACCTACTTCGGGCCATCCTGAATTCCTCAGGC
 ATGGAAAATTGGTGTGACAAGGACATGGCCGAGGAGGGGTCCTGGAGGAAGCAGAGTTCTAACACAT
 TGGTCTCTGATTGCAATCATCAAAGACAGGATGGAGGAGAAAGACTACACAGTGGCCAGGTGCCACCT
 AAGCATGTATACCGAGTGTGACGTGCCAAGAGGAGGAGCTCACACAGATGGTATCCACTATGTCCGACG
 GCTGGCGCTTCGAGCAGCTGGTGAACATCGGCTCCTCCTACAACATATGGGAGTGAGGACCAGGCCGAGTT
 CCTGTGTGTGGTGTCTAAGGAGCTCCACAGCTCCCCACACGGCCTGAGCTCCGAGTCCACTCGCAAGGCC
 AAGCCCTGTATCTCAGGATCCCGCTAACCTTTCTCCCTCCCACCACCGCCTCCTCCTCCGCTTCC
 TGCTGGAGTCCCTGTCCGACCCTCCAGACCCAAGCCTGAGCTTGAGTGA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001289671
Insert Size: 753 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001289671.1, NP_001276600.1</u>
RefSeq Size:	1707 bp
RefSeq ORF:	753 bp
Locus ID:	72844
Cytogenetics:	15 E1