

Product datasheet for **MC226612**

Tnfrsf9 (NM_011612) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Tnfrsf9 (NM_011612) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Tnfrsf9
Synonyms:	4-1BB; A930040I11Rik; AA408498; AI325004; Cd137; CDw137; ILA; Ly63
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>MC226612 representing NM_011612 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGGAAACAACGTGTACAACGTGGTGGTTCATTGTGCTGCTGCTAGTGGGCTGTGAGAAGGTGGGAGCCG
TGCAGAACTCCTGTGATAAAGTGCAGCCTGGTACTTTCTGCAGAAAATACAATCCAGTCTGCAAGAGCTG
CCCTCAAAGTACCTTCTCCAGCATAGGTGGACAGCCGAAGTGAACATCTGCAGAGTGTGTGAGGCTAT
TTCAGGTTCAAGAAGTTTGTCTCTACCCACAACGCGGAGTGTGAGTGCATTGAAGGATTCATTGCT
TGGGGCCACAGTGCACCAGATGTGAAAAGGACTGCAGGCCTGGCCAGGAGCTAACGAAGCAGGGTTGCAA
AACCTGTAGCTTGGGAACATTTAATGACCAGAACGGTACTGGCGTCTGTGACCCCTGGACGAACTGCTCT
CTAGACGGAAGGTCTGTGCTTAAGACCGGGACACGGAGAAGGACGTGGTGTGTGGACCCCTGTGGTGA
GCTTCTCTCCAGTACCACCATTTCTGTGACTCCAGAGGGAGGACCAGGAGGGCACTCCTTGAGGTCCT
TACCTTGTTCCTGGCGCTGACATCGGCTTTGCTGCTGGCCCTGATCTTCACTCTCCTGTTCTCTGTG
CTCAAATGGATCAGGAAAAAATCCCCACATATTCAAGCAACATTTAAGAAGACCACTGGAGCAGCTC
AAGAGGAAGATGCTTGTAGCTGCCGATGTCCACAGGAAGAAGAAGGAGGAGGAGGCTATGAGCTGTG
A

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	Sgfl-MluI
ACCN:	NM_011612
Insert Size:	771 bp



[View online »](#)

OTI Disclaimer:	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info</p>
OTI Annotation:	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_011612.2, NP_035742.1</u>
RefSeq Size:	2134 bp
RefSeq ORF:	771 bp
Locus ID:	21942
UniProt ID:	<u>P20334</u>
Cytogenetics:	4 81.52 cM
Gene Summary:	<p>Receptor for TNFSF9/4-1BBL. Possibly active during T cell activation.[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (1) represents the longest transcript. Variants 1 and 3 encode the same protein (isoform 1).</p>