

## Product datasheet for **MC226648**

### Vps37d (NM\_001199677) Mouse Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Vps37d (NM\_001199677) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Vps37d  
**Synonyms:** Wbscr24  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >MC226648 representing NM\_001199677  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**GCGATCGC**C

ATGTACCGGGCCCGGGCGGGCGGGCGGGCGGGAGCCCGGAGCCCGGGCGCTTTGGGATCCTCAGCA  
CCGGGCAGCTCCGGGATCTGCTTCAAGATGAGCCCAAGCTGGACCGGATCGTGCGGCTCAGCAGGAAGTT  
CCAGGGTCTGCAGCTGGAACGGGATGCCTGCTTGGCCTCTAACTACGCGCTGGCCAAGGAGAACCCTGGCC  
TTGCGGCCCGCCTGGAGATGGGTCGGACAGCCCTGGCCATCAAGTACCAAGAGCTTCGAGAGGTGGCTG  
AGAAGTGCAGACAACTGCAGAGACTGGAGAAGAGCATGCATCGCTGGAGCCCAATGTGCGCTGGG  
CTGGCTGCAGGCTGAGCTGGAAGAGGCTGAGCAGGAGGCTGAGGTCCAGATGGAGCAGCTGCTGCTGGGG  
GAACAGAGCCTGGAGGCCTTCTGCCAGCCTTCCAGAGGGCCGTGCCCTGGCCACCTGAGGCGGACAC  
AGGCCGAGAACTACAGGAGGTGCTGAGGCGTCGAGAGCGCTCTGCCAGCCGGCCCCACAACCTGCTGC  
TGCTGCTGCCCGCCGCCACCGCCATGGATCCCCCAAACCTTCCCTGCTGCTGCCGCTCTGCCACT  
GGGGCCCGCCGGGACCACCTCCAGCAGTACCCCGGAGCCTGCCCCCTTGGACTCCCGCCAGTGCCCC  
CAGTGAAGGGCTCTCCGGGTGCCCTTTGGCCAGCCCTTTGCTGAGCCCTCGACCCTCAGCCTGA  
ACCCCTCATCG**TAG**

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_001199677  
**Insert Size:** 786 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u>NM_001199677.1, NP_001186606.1</u>
<b>RefSeq Size:</b>	1561 bp
<b>RefSeq ORF:</b>	786 bp
<b>Locus ID:</b>	194309
<b>UniProt ID:</b>	<u>Q810I0</u>
<b>Cytogenetics:</b>	5 G2
<b>Gene Summary:</b>	Component of the ESCRT-I complex, a regulator of vesicular trafficking process. Required for the sorting of endocytic ubiquitinated cargos into multivesicular bodies. May be involved in cell growth and differentiation (By similarity).[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (1) represents the longer transcript and it encodes the longer protein (isoform 1).