

## Product datasheet for **MC226666**

### Gm20765 (NM\_001270644) Mouse Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Gm20765 (NM\_001270644) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Gm20765  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >MC226666 representing NM\_001270644  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGACTACCATGGCTGACAAAAATGGACATGTCATTGGAAGACATCATTAAAGCTGACCAAGATCCAGCAAA  
GACGTCATGACAGGCCTGATCCCGGGTGAATCGAGGCACTGGGTCTAAGAGGTATAGGCCAGCTTTCAC  
CCATGGTGAAGGAACAGGCTGGCACCATACTGCAGACCAAAACAACCTCCAGACAAGTGGCAGCATGAC  
CTGTTTATTGGAGGCTTCAGGGGTCAAAACCACGTGGATACTGGTGGGAAGCTCTTCTTGTCAAATCTGC  
ACTTCGGAGTGTGAGATGCTGATATTGAGTACTCTTTGCTGAATTTGGAACGTTGAAGAAATCTGCTGT  
GCACTATGATCGCTGTGGACGAAGTTTAGGGACAGCACAGGTGCACCTTTGAAAGGAAAGCAGATGCCCTG  
AAGGCTATGAGAGAGTACAATGGCGCCCTTTGGATGGCCGCCCTATGAACATCCAGCTTGCCACCTCAC  
AGATTGATAGACAAGGAAGACCTGCACAAAGCAAAAATAGGGGCGGCATGACAAGAAACCCTGGCTCTGG  
AGTATTAAGTGGTGGAGGCACCAAGAAATGGACACTTGGAGGCAGCCAGGGAAGAGGGAGAGGCACCATC  
AGGAACTCAAACTGCAGCAACAGCAGCAGCAGCAGAGAAGCAGAAGCAGCAGCAGCAGCAGAGAAGCAGC  
AGCAGAAGCAGCAGCAGAAGCAGCAGCAGCAGCTATCTGCAGAGGAGCTGGATGCCAGCTGGATGCTTA  
TCAGGAAATGATGGACACCAGCTGA

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_001270644  
**Insert Size:** 795 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001270644.1</a></u> , <u><a href="#">NP_001257573.1</a></u>
<b>RefSeq Size:</b>	1246 bp
<b>RefSeq ORF:</b>	795 bp
<b>Locus ID:</b>	637578
<b>Cytogenetics:</b>	10 D1