

Product datasheet for **MC226695**

Ms4a7 (NM_001276398) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Ms4a7 (NM_001276398) Mouse Untagged Clone
Tag: Tag Free
Symbol: Ms4a7
Synonyms: 9130422I10Rik; A430103C15Rik; CD20I4; CFFMA
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC226695 representing NM_001276398
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGCTTCAGATACAACCAACATTTGGCTTCCAATGCAGTGCCATCCCCATCTGCCCCACTGGTAAACAAGA
ATCCGGGGTTCATTTGCATCCCTGGCACCATGCGTCTACAGCTTGGCACCAAGAACATTGGGTGGGACTG
CTTTCCAAAGGACATCATTATCCACAAAAGAGAGAAAACCTGGACATACATATGAAAAAGAAGATGACCTG
CTGATTGGAGTGCCTAGTGAAGCCACAGTTCTTGGAAACATCCAGCTACTATGTGCCCTGATACTTGCAA
GCTTTGGGGGCATTTTGGTGTGTCAGCTTCTACTTCAACCCAGAAGTTCCACCACCTTTGATATCTGGATA
CCTGTTTATTGGATCCCTCTGTTTTGCCATTGCTGGAATTCTCTCCATTATCTCTGAAAAATCTCAACC
AAACCTTTGCACTGAGCAGCCTGGCCTCCAATGTAGCAAGCTCTGTTGTTGCCGTATTGGCCTCTTCC
TCTTCACTATTGTCTGATAGCCCTGGGGAGTGCCTTCCACACTGTAACCTCAGAAAAGAAATCCTGTC
CTTATTGTCTTACTTGAAATCCCACCATTTGGAAGAATGAAGACAAGAAGTGTACCTGGCTTATGTCGGT
GCTATGAGTGTCTGGGGATGATGCTTCTTCACTGTGCTGGAGGTTTCTAGCTGGATACAGTTCTA
TCTTTTGGTGAAGCAAGTCTACTCCAACAACTGGGGGTACATTTTCTTGCCTCAATCACAAGATCA
TACCCAACCTGTTAAAAGCAATTTGTTGCAATAA

ACGCGTACGCGGCCGCTCGAGCAGAAAACCTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001276398
Insert Size: 804 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001276398.1</u> , <u>NP_001263327.1</u>
RefSeq Size:	1576 bp
RefSeq ORF:	804 bp
Locus ID:	109225
Cytogenetics:	19 A