

## Product datasheet for **MC226710**

### Ube2j2 (NM\_001284312) Mouse Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Ube2j2 (NM\_001284312) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Ube2j2  
**Synonyms:** 1200007B18Rik; 2400008G19Rik; 5730472G04Rik; AL022923; NCUBE-2; Ubc6; Ubc6p  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >MC226710 representing NM\_001284312  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**GCGATCGCC**

ATGCATCAGGCCAAAAGCCCTTGATTCTGCCAGAAAGATGAGCAATAACAGCAATAAGAGAGCTCCAACAA  
CAGCAACACAAAAGTTGAAGCAGGACTACCTTCGAATTAAGAAAGACCCAGTGCCTTACATCTGTGCAGA  
GCCCTCCCTCAAATATTCTTGAATGGCATTATGTTGTCCGAGGCCCTGAGATGACTCCTTATGAAGGT  
GGCTATTATCATGGAAAATAATTTTTCCAGAGAATTTCCATTTAAGCCTCCTAGTATTTACATGATAA  
CCCCAATGGAAGATTTAAGTGCAACACAAGGCTGTGTCTCTCCATCACAGATTTCCACCCAGATACATG  
GAATCCAGCGTGGTCCGCTCCACCATTCTAACAGGGCTCCTGAGCTTTATGGTGGAGAAAGGCCCCACC  
CTGGGGAGCATAGAGACGTCGGACTTCACGAAAAACAGCTGGCAGCACAGAGCTAGTGTTTAATTTAA  
AAGATAAAGTCTTTTGTGAATTTTCTGAAGTTGTGGAGGAAATTAACAGAAACAGAAAGCACAAGA  
TGAACAAAGTAAACAGGCCCCAGAATCTTCTTTGCCAGACGTGGTTCCTGATGGGGAGCTGCACCGTGGC  
CAGCATGGGATCCAGCTCCTCAATGGGCACGCACCCGCTGCTGGCCCCAACCTCGCTGGGCTCCCACAGG  
CCAACAGGCATCACGGACTCCTGGCGGGCGCTTGGCGAACTTGTGGTATAGTTGGGTTTGCAGCCTT  
CGCCTACACGGTCAAGTACGTGCTGAGAAGCATAGCACAGGAGTGA

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_001284312  
**Insert Size:** 816 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u>NM_001284312.1, NP_001271241.1</u>
<b>RefSeq Size:</b>	3284 bp
<b>RefSeq ORF:</b>	816 bp
<b>Locus ID:</b>	140499
<b>UniProt ID:</b>	<u>Q6P073</u>
<b>Cytogenetics:</b>	4 87.66 cM
<b>Gene Summary:</b>	<p>Catalyzes the covalent attachment of ubiquitin to other proteins. Seems to function in the selective degradation of misfolded membrane proteins from the endoplasmic reticulum (ERAD).[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (5) differs in the 5' UTR compared to variant 1. Variants 1, 2 and 5 encode the same isoform (a). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>