

Product datasheet for MC226770

Sult2a8 (NM_001199306) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids

Tag: Tag Free

Symbol: Sult2a8

Synonyms: 2810007J24Rik; mL-STL

Mammalian Cell Selection: Neomycin

Vector: pCMV6-Entry (PS100001)

E. coli Selection: Kanamycin (25 ug/mL)

Fully Sequenced ORF: >MC226770 representing NM_001199306
Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTGAATACGACTCACTATAGGGCGGCCGGAATTCGTGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGACAGATGAATTTCTGTGGATAGAAGGGATACCATTCCCTACAGTTTACTACTCACAAGAAATCATAA
GAGAAGTACGTGACAGGTTTGTGGTAAGGGATGAAGACACGATCATAGTGACTTATCCAAAATCAGGAAC
CCACTGGTTGAATGAAATTGTCTGCTTGATTCTGACCAAGGGAGATCCACATGGGTCCAATCTACAATT
GCCAATGAGCGCACACCATGGATAGAATTTGAAAATAACTATAGAATACTAAACAGCAAGGAGGGTCCAC
GTCTCATGGCCTCTCTCCTCCCATCCAATTTCCAAAGTCTTTCTTCAGTTCTAAGGCAAAGGTGAT
TTATCTCATCAGAAATCCAGAGATGTTCTTGTGTGAGGCTATCACTACTTCAATGCTCTAAACAGGGA
AAGGAACAAGTGCCATGAAAAATATTTTTGAGAACTTTTGAAGGAAATCATATTTGGATCATGGT
TCGAACATGCATGTGGCTGGATATCCCTGAGGAAAAGAGAAAACATTCTAGTGTTAAGTTATGAACAGCT
GAAAAAGGATACAAGAAATACCATCAAGAAGATCTGTGAATTCCTGGGAGAGAAATTTGGAATCAGGGGAA
CTGGAGTTAGTTTTGAAGAACATTTCTTTTCAGATCATGAAAGAAAGGATGATATCTCAGTCATGTCTCT
CAAATATAGAAAAACATGAATTCATTATGAGAAAAGGCATCACAGGGGACTGGAAGAATCACTTTACAGT
AGCCCAGGCTGAAGCCTTTGATAAAGCATTTTCAGGAGAAAGCAGCTGATTTTCTCAAGAGCTCTTTTCA
TGGGAATAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI

ACCN: NM_001199306



Insert Size:	849 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	NM_001199306.1 , NP_001186235.1
RefSeq Size:	2097 bp
RefSeq ORF:	849 bp
Locus ID:	76971
Cytogenetics:	7 A1