

Product datasheet for **MC226811**

Otx2 (NM_001286482) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Otx2 (NM_001286482) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Otx2
Synonyms:	E130306E05Rik
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>MC226811 representing NM_001286482 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**

ATGATGTCTTATCTAAAGCAACCGCCTTACGCAGTCAATGGGCTGAGTCTGACCACTTCGGGTATGGACT
TGCTGCATCCCTCCGTGGGCTACCCGCCACCCCGGAAACAGCGAAGGGAGAGGACGACATTTACTAG
GGCACAGCTCGAGTTCTGGAAGCTCTGTTTGCCAAGACCCGGTACCCAGACATCTTCATGAGGGAAGAG
GTGGCACTGAAAATCAACTTGCCAGAATCCAGGGTGCAGGTATGGTTTAAAGATCGAAGAGCTAAGTGCC
GCCAACAGCAGCAGCAGCAGCAGAATGGAGGTGAGAACAAGTGAGGCCTGCCAAGAAGAAGAGCTCTCC
AGCTCGGGAAGTGAGTTCAGAGAGTGGAACAAGTGCCAGTTCAGTCCCCCTCTAGTACCTCAGTCCCA
ACCATTGCCAGCAGCAGTCTCAGTGTCTATCTGGAGCCAGCGTCCATCTCCCACTGTCTGACCCCT
TGTCCACTTCCTCCTCCTGCATGCAGAGGTCTATCCCATGACCTATACTCAGGCTTCAGGTTATAGTCA
AGGCTATGCTGGCTCAACTTCCTACTTTGGGGCATGGACTGTGGATCTTATTTGACCCCTATGCATCAC
CAGCTTCCTGGACCAGGGGCCACACTCAGTCCCATGGGTACCAATGCTGTTACCAGCCATCTCAATCAGT
CCCCAGCTTCTTTCCACCCAGGGATATGGAGCTTCAAGCTTGGGTTTTAACTCAACCACTGATTGCTT
GGATTATAAGGACCAAACTGCCTCTTGAAGCTTAACTTCAATGCTGACTGCTTGGATTATAAGATCAG
ACGTCTCATGGAAATTCAGGTTTT**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	SgfI-MluI
ACCN:	NM_001286482
Insert Size:	870 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_001286482.1</u> , <u>NP_001273411.1</u>
RefSeq Size:	2317 bp
RefSeq ORF:	870 bp
Locus ID:	18424
Cytogenetics:	14 25.36 cM
Gene Summary:	This gene encodes a protein that belongs to the homeobox family of transcription factors. The encoded protein plays a role in the development and patterning of the head. This protein regulates development of the choroid plexuses in the brain affecting composition of cerebrospinal fluid in the developing brain and is thought to function in the development of sense organs in the embryo. In humans, mutations in this gene are associated with pituitary hormone deficiency. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Nov 2013] Transcript Variant: This variant (2) uses an alternate in-frame splice site in the 5' coding region compared to variant 1. The encoded protein (isoform b) is shorter than isoform a. Variants 2, 3 and 4 encode the same protein.