

## Product datasheet for **MC226830**

### Kansl2 (NM\_001289439) Mouse Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Kansl2 (NM\_001289439) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Kansl2  
**Synonyms:** 2310037124Rik; B930009D17  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >MC226830 representing NM\_001289439  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**GCGATCGCC**

ATGAGAGAAAAGCTAATTCGTTTGCAGTCTTTGTACATTGATCAGTTTAAAAGACTCCAACACCTGCTCA  
AAGAGAAGAAGCGCCGCTATTTACACAACCGGAAAGTGGAGCACGAGGCTTTAGGTAGTAGTCTCCTGAC  
TGGCCCTGAGGGACTTTTGGCCAAAGAGAGAGAGAATTTAAAGCGACTAAAATGTCTTAGCGGATACCGT  
CAACGCTATGGAGTAGAAGCCTTACTACACAGGCAGCTGAAGGAACGCAGAAATGCTGGCCACCGATGGGG  
CTGCTCAACAGGCTCACACCACTCGTTCAGTCAGAGGTGCTTAGCCTTTGTGGATGATGTGCGCTGCTC  
CAATCAGTCTCTCCCAATGACCAGACACTGCCTTACCCATATCTGTCCAGGATACAATCAGGTTCTCTTC  
AAATGCTGCCAGGGATCTGAAGAGGTCCCTGCAACAAACAGTTCTCTGTAAGCCTTTCTGAGGATCCCT  
GCTGCCCTCTGCATTTCCAGTTGCCTCCTCAGATGTATAAGCCCGAGCAGGACTGTCTGTGCCAGACGA  
TCTGGAAGCCGGCCCATGGATCTGTACTTGAGTGCTGCCGAGCTTCAGCCCACTGAGAGTTTACCAGT  
GAGTTCAGTGATGACTTGGATGTTGTGGGTGATGGCATGCCATGCCCGCCGTCGCCATTGCTCTTTGATC  
CTTCGCTAACTCTTGAAGATCATTAGTACAGAAAGTGCCGAGACGCAGGACAGATGCAGGCGGCTGG  
GGATGGGTGCAGATCCCAGGGAGCTCGAAGTGTAGAGAAAGCAGCATTCCTCAGCGTGGATTGGCTACT  
GCGAATGGAAAACCAAGAACCACTTCTATCAGT**TGA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_001289439  
**Insert Size:** 876 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001289439.1</a></u> , <u><a href="#">NP_001276368.1</a></u>
<b>RefSeq Size:</b>	3502 bp
<b>RefSeq ORF:</b>	876 bp
<b>Locus ID:</b>	69612
<b>UniProt ID:</b>	<u><a href="#">Q8BQR4</a></u>
<b>Cytogenetics:</b>	15 F1
<b>Gene Summary:</b>	<p>As part of the NSL complex it is involved in acetylation of nucleosomal histone H4 on several lysine residues and therefore may be involved in the regulation of transcription.</p> <p>[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (4) lacks an alternate exon compared to variant 1. The resulting isoform (d) is shorter at the N-terminus compared to isoform a. Variants 4 and 5 both encode the same isoform (d).</p>