

## Product datasheet for **MC226851**

### Nipa2 (NM\_001256133) Mouse Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Nipa2 (NM\_001256133) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Nipa2  
**Synonyms:** 2600017P10Rik; 3830408P04Rik; AB041581  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >MC226851 representing NM\_001256133  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

ATGGGAGCCGGTGAAGTGGCCAATTTTGTGCATATGCATTTGCTCCGGCCACATTAGTGACTCCACTAG  
GTGCGCTCAGCGTCTCGTAAGTGCCATTCTGTCTTCATACTTCTAAATGAAAGACTTAATCTTCATGG  
GAAAATTGGGTGTTTGTAAAGTATTCTAGGGTCTACAGTTATGGTCATTCATGCTCCAAAGGAAGAGGAG  
ATTGAGACCTTAAATGAAATGTCACATAAGTTAGGTGATCCAGGTTTTGTGGTCTTTGCAACATTTGTGG  
TCATTGTGGCCTTGATATTCATCTTTGTGGTGGGACCTCGACATGGACAGACAAACATTTCTGTGTATAT  
AACAACTGCTCTGTGATTGGAGCATTTCAGTCTCCTGTGTGAAGGGTTTGGGTATCGCCATAAAAGAG  
TTGTTGGCTGGAAAGCCTGTTCTGCAGCACCCCTCGCCTGGATTCTTCTGTTTAGCCTTGTGGTCTGTG  
TGAGCACACAGATTAATTACCTAAACCGGGCTCTGGACATTTTCAACACTTCCATTGTGACTCCAATATA  
TTATGTGTTCTTTACAACCTCAGTTTTGACTTGTTCAGCTATTCTTTTAAAGAATGGCAGGATATGCCT  
GTTGATGATGTCATTGGGACTTTGAGTGGCTTCTTCAACATCATTGTGGGATATTTTGTGTCATGCCT  
TTAAAGATGTCAGTTTTAGTCTAGCAAGTCTGCCAGTATCTTTTCGAAAAGATGAAAAGCAATGAATGG  
CAATCTTTCTAGTATGTATGAAGTCTTAATAATAATGAAGACGACTTACCCTGTGGAATTGAACACACT  
GGTGAAAATATTTCCCGTAGGAATGAAAATCTGCCATCTTTTT**TAA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_001256133  
**Insert Size:** 885 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001256133.1</a></u> , <u><a href="#">NP_001243062.1</a></u>
<b>RefSeq Size:</b>	2723 bp
<b>RefSeq ORF:</b>	885 bp
<b>Locus ID:</b>	93790
<b>UniProt ID:</b>	<u><a href="#">Q9JJC8</a></u>
<b>Cytogenetics:</b>	7 B5
<b>Gene Summary:</b>	Acts as a selective Mg(2+) transporter.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (5) differs in the 5' UTR, lacks a portion of the 5' coding region, and initiates translation at a downstream AUG start codon, compared to variant 1. The encoded isoform (2) has a shorter N-terminus, compared to isoform 1.