

Product datasheet for MC226866

Sirpa (NM_001291021) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Sirpa (NM_001291021) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Sirpa
Synonyms:	AI835480; Bit; CD172a; P84; Ptpns1; SHP-1; SHPS-1; SIRP
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>MC226866 representing NM_001291021 Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTGAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCCGCATCGCC

ATGGAGCCCGCCGCGCCCGCCCTGGCCGCTAGGGCCGCTGCTGCTCTGCCTGCTGCTCTCCGCGTCT
GTTTCTGTACAGGAGCCACGGGAAGGAAGTGAAGGTGACTCAGCCTGAGAAATCAGTGTCTGTTGCTGC
TGGGGATTCGACCGTTCTGAACGCACTTTGACCTCCTGTTGCCGTGGGACCCATTAGGTGGTACAGA
GGAGTAGGGCCAAGCCGCTGTTGATCTACAGTTTCGAGGAGAATACGTTCTCGAATTAGAAATGTTT
CAGATACTACTAAGAGAAACAATATGGACTTTTCCATCCGTATCAGTAATGTACCCAGCAGATGCTGG
CATCTACTACTGTGTGAAGTTCCAGAAAGGATCATCAGAGCCTGACACAGAAATACAATCTGGAGGGGA
ACAGAGGTCTATGTAATCGATAAATGTACCCACAAGTGAATGTCTTCATCGGTGTGGCGTGGCGT
GTGCTTTGCTCGTAGTCTGCTGATGGCTGCTCTACCTCCTCCGGATCAAACAGAAAGGCCAAGGG
GTCAACATCTTCCACACGGTTGCACGAGCCGAGAAGACGCCAGGGAAATAACCCAGGTACAGTCTTTG
ATCCAGGACACAAATGACATCAACGACATCACATACGACAGCTGAATCTGCCAAAGAGAAGAGCCCG
CACCCCGGGCCCTGAGCCTAACACCACAGAAATATGCAAGCATTGAGACAGGCAAGTGCCTAGGCC
AGAGGATACCTCACCTATGCTGACCTGGACATGGTCCACCTCAGCCGGGCACAGCCAGCCCCAAGCCT
GAGCCATCTTTCTCAGAGTATGCTAGTGTCCAGGTCCAGAGGAAGTGA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	SgfI-MluI
ACCN:	NM_001291021
Insert Size:	888 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_001291021.1</u> , <u>NP_001277950.1</u>
RefSeq Size:	3389 bp
RefSeq ORF:	888 bp
Locus ID:	19261
UniProt ID:	<u>P97797</u>
Cytogenetics:	2 63.19 cM
Gene Summary:	<p>Immunoglobulin-like cell surface receptor for CD47. Acts as docking protein and induces translocation of PTPN6, PTPN11 and other binding partners from the cytosol to the plasma membrane. Supports adhesion of cerebellar neurons, neurite outgrowth and glial cell attachment. May play a key role in intracellular signaling during synaptogenesis and in synaptic function. Involved in the negative regulation of receptor tyrosine kinase-coupled cellular responses induced by cell adhesion, growth factors or insulin. Mediates negative regulation of phagocytosis, mast cell activation and dendritic cell activation. CD47 binding prevents maturation of immature dendritic cells and inhibits cytokine production by mature dendritic cells (By similarity).[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (6) lacks two alternate in-frame exons in the central coding region, compared to variant 4, resulting in an isoform (4) that is shorter than isoform 3.</p>