

## Product datasheet for **MC226871**

### Fam110a (NM\_001289150) Mouse Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Fam110a (NM_001289150) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Fam110a
Synonyms:	1700008J10Rik; 5430432M24Rik; AI747080
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
Fully Sequenced ORF:	>MC226871 representing NM_001289150 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**GCGATCGCC**

ATGCCTGTGGATACGCTGAGCCCCGGAGCCCCGCCACTCCCGCCCTACCTTTTCGCTGAGAACCAAGG  
TTCCCGGCTACCTGCTACCAAGGCCAGCAGACGGGGGAGCCAGGAAACCGAGTGTGTGGAACGCCTGGA  
GGCTGATAAGGCTAAGTACGTCAAGAGCCTGCGTGTGGCCAACACCCGCCAAGAGCCGGTGCAGCCCCCG  
CTGGTCAGACAGCCCCTCTCAGCCCTGGTCTCGCGGCCAGTGTCTCACACCCAGTCGCCGAGTCTCTGC  
CCTGTTCTGGCCGCCGGCCCCAACTGGACCTCGACATCCTTAGCAGCCTCATCAACTTGTGTGATAGTCC  
TGTGTCCCATCTGAGGCCAGCCGAACACTGGAAGGCCAGAAGGATCTGCCACAAGTCCCACCAGCC  
ACCCCTCCTCGTCCGCCACCCAGTACTGTGCGGGTGGCCGAGTGGACGTTTCGTCCTTGCCCTGCTTAC  
CTGCTCGGCCCTATCCATCACCTGGCACCCTACCCTTCCAGCCCCGGCCGGCCCCCTGGCCTGCAGCG  
GTCCAAGTCGGACTTGAGTGAGCGCTTTCCAGGGCAGCAGCCGATCTTGAGCGCTTTTTAACTTCTGT  
GGCCTGGACCCAGAGGAAGCGCGAGGATTGGGGTGGCCACTTAGCAAGGGCCAGCTCGGACATCGTGT  
CTCTAGCCGGCCAAGTGTGGACCTTGACGCTGAAGGGGGCTGCTCACGCCGAGCTCTGCTACAGT  
GGAAGAGCGGTCCCTGGATCGTGTCCCCTATGGGGTGTCTGTGATCGAGCGAAACGCTCGTGTGATCAAG  
TGGCTATATGGGTTGCGGCAGGCTCGTGACCCTCCCACCACTGAGGGT**TAG**

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	Sgfl-MluI
ACCN:	NM_001289150
Insert Size:	891 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001289150.1</a></u> , <u><a href="#">NP_001276079.1</a></u>
<b>RefSeq Size:</b>	1608 bp
<b>RefSeq ORF:</b>	891 bp
<b>Locus ID:</b>	73847
<b>UniProt ID:</b>	<u><a href="#">Q8R184</a></u>
<b>Cytogenetics:</b>	2 G3