

## Product datasheet for **MC226880**

### Tmem241 (NM\_001289666) Mouse Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Tmem241 (NM_001289666) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Tmem241
Synonyms:	6030446N20Rik
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>MC226880 representing NM_001289666 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGAACGTGAGGCGGTCCCTGCTCGGCCTCACCTTCTGCACCTGTTACCTGGTCTCACCTACCAACA  
AGTATGTCCTGTCTGCTGAAATTTACCTACCCTACATTATCCAAGGGTGGCAGACATTCATTGGTGG  
ACTTTTGTTCATATGTCATGGAACTAGGATGGGTGGAGCTCCACAGCAGTCCAAGATCTGATGTTTTG  
ATATGGCTTCTGCTTCAGCGCTGTTTGTGGGCATAATCTATGCGGGTCCAAGCTTTGTCTAGACTGG  
CTGTTCTGTGTTTTTCATTTTGCATAATGTAGCTGAAGTCTCACCTGTGGGTACCAGAAGTGTGTCTG  
GAAAGAGAAAACATCGCTCTCAAAGATCTGTAGTGCCCTCTTCTCTTGGCTGCTGCAGGGTGCCTTCCC  
TTCCAAGACTCTCAGTTTGATCCAGATGGATATTTTTGGGCTCTTATCCACATCTTTTGTGTAGGTTCTT  
ACAAAATACTTCGGAAGTCCCGGAAGCCCACTGTGCTAAGTGACATTGACCAGCAGTACTTAAACTATAT  
ATTCAGCATGGTACTCCTGGCATTTCATCTCATCCGACAGGTGACCTCTTCGGTGCCTGGACTTTCCA  
TTCTTGTATTTCTACAGGTTCCATGGCAGCTGCTGTGCCAGTGGAGTTTTGGGATTTCTCCTCATGCTCA  
GTACAGTCAGGCTGAGAAGCATCCTGGCCCCAGGGCAGTGGCAGCCTGGATCCTCTGTGCTAAGGTCGT  
CACAGCTGGCTTGTGATGTTGCTGTTTGTATATGGCCCTGACCAAAGCAACTGTGGGATGCTTCTGCTT  
GGTGGACTTGGCAGGCCCTGCTGGTTTTCTCAGAGCGGAGGAGTTCCTCA**TGA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	Sgfl-MluI
ACCN:	NM_001289666
Insert Size:	894 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001289666.1</a></u> , <u><a href="#">NP_001276595.1</a></u>
<b>RefSeq Size:</b>	3206 bp
<b>RefSeq ORF:</b>	894 bp
<b>Locus ID:</b>	338363
<b>UniProt ID:</b>	<u><a href="#">Q3UME2</a></u>
<b>Cytogenetics:</b>	18 A1