

## Product datasheet for **MC226901**

### Sh3yl1 (NM\_001289480) Mouse Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Sh3yl1 (NM\_001289480) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Sh3yl1  
**Synonyms:** AI314953; Ray; YSC84  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >MC226901 representing NM\_001289480  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**GCGATCGCC**

ATGAATAACCCGATACCTTCTAATCTGAAATCAGAAGCAAAGAAGGCTGCCAAAATATTAAGAGAATTCA  
CAGAAATAACTTCTAGAAATGGACCTGATAAGATCATTCTGAATGGTCTGCACCTTCAGCCATTGGGAT  
AGCTGGGCTTGGTGGAGGATTTGAAATAGGAATAGAGGTATCCGATTTGGTGATAATTCTGAATTATGAC  
AGAGCAGTAGAAGCCTTTGCAAAGGGTGGAAACCTGACCCTGGAGGGAATTTACTGTCCGAGTTGGGC  
CCCTAGGCAGGAACCTGGAAGGCAACGTGTCCCTGAGGAGCTCCGCAGCCGTCTTACATACTGCAAGTC  
CAGGGGCTCTTTGCTGGTATCTCCCTGGAAGGAAGCTGTCTGATTGAGAGGAAAGAAACCAATCGGAAA  
TTTTACTGTCAGGACATTCGAGCTTACGACATTTATTTGGAGATGTTCCGCAGCCTGCTCAAGCAGAAG  
ATCTTTATGAAATTTAATTCCTTCACTGAGAAATATGAACTGAAGGGCAGCGCATCAATTTGAAGAA  
AGTAGCCCGGGAGCAGAGGAAGGCCAAAGAATTACCTCCAAAACCATCTCAAGGCCCCAGCCAGCCAC  
CCACCTGTACAGCTGAATGCTGGCTCCCAGGTAATAGAAATGAATATAAGCTATATCTGAGCTTTCTA  
GCTATCATGAGAAACTGGCAATTTGAATCAACCAATAGAAGTGACAGCATTATACTCTTTCCGAAAGGACA  
ACAGCCAGGAGACTTGAATTTTCAAGCAGGAGACAGAATCATAGTTATATCAAAGACCGACTCAAATTTT  
GATTGGTGGGAAGGAAAACCTTCGAGGTCAAACCTGGAATTTCCAGCCAACCTATGTAACCTATGAAT**TAA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAAAACCTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_001289480  
**Insert Size:** 909 bp



[View online »](#)

|                               |   |
|-------------------------------|---|
| <b>OTI Disclaimer:</b>        | Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).  |
| <b>OTI Annotation:</b>        | Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.  |
| <b>Components:</b>            | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).  |
| <b>Reconstitution Method:</b> | <ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol> |
| <b>RefSeq:</b>                | <u><a href="#">NM_001289480.1</a></u> , <u><a href="#">NP_001276409.1</a></u>   |
| <b>RefSeq Size:</b>           | 1633 bp   |
| <b>RefSeq ORF:</b>            | 909 bp  |
| <b>Locus ID:</b>              | 24057   |
| <b>UniProt ID:</b>            | <u><a href="#">O08641</a></u>   |
| <b>Cytogenetics:</b>          | 12 A2   |