

## Product datasheet for **MC226949**

### Ager (NM\_001271424) Mouse Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Ager (NM\_001271424) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Ager  
**Synonyms:** RAGE  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >MC226949 representing NM\_001271424  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGC**C

ATGGAAGTCCACTGGAATTGTCGATGAGGGGACTTCCGGTGTGGGCAACTAACAGGCGAGGGAAGG  
AGGTCAAGTCCAACCTACCGAGTCCGAGTCTACCAGATTCTGGGAAGCCAGAAATGTGGATCCTGCCTC  
TGAAGTCAAGCCAGTGTCCCTAATAAGGTGGGACATGTGTGTCTGAGGGAAGCTACCTGCAGGGACC  
CTTAGCTGGCACTTAGATGGGAACTTCTGATTCCCGATGGCAAAGAACTCGTGAAGGAAGAGACCA  
GGAGACACCCTGAGACGGGACTTTTACTGCGGTGACAGCTGACAGTATCCCCACCCAAAGGAGGAAC  
CCATCCTACCTTCTCTGCAATTCAGCCTGGGCTTCCCGGCGCAGACCCCTGAACACAGCCCCATC  
CAACTCCGAGTCAGGGAGCCTGGGCTCCAGAGGGCATTAGCTGTTGGTTGAGCCTGAAGGTGGAATAG  
TCGCTCCTGGTGGGACTGTGACCTTGACCTGTGCCATCTGCCCAGCCCCCTCCTCAGGTCCACTGGAT  
AAAGGATGGTGCACCCTTGCCCCGGCTCCCAGCCCTGTGCTGCTCCTCCCTGAGGTGGGGCAGCAGGAT  
GAGGGCACCTATAGCTGCGTGGCCACCCACCCTAGCCACGGACCTCAGGAAAGCCCTCCTGTCAGCATCA  
GGTTCACAGGCTCTGTGGGTGAGTCTGGGCTGGGTACGCTAGCCCTGGCCTTGGGGATCCTGGGAGCCT  
GGGAGTAGTAGCCCTGCTCGTGGGGCTATCCTGTGGCGAAAACGACAACCCAGGCGTGAGGAGAGGAAG  
GCCCGGAAAGCCAGGAGGATGAGGAGGAACGTGACAGCTGAATCAGTCAGAGGAAGCGGAGATGCCAG  
AGAATGGTCCCGGGGACCG**TAA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_001271424  
**Insert Size:** 933 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001271424.1</a></u> , <u><a href="#">NP_001258353.1</a></u>
<b>RefSeq Size:</b>	1455 bp
<b>RefSeq ORF:</b>	933 bp
<b>Locus ID:</b>	11596
<b>Cytogenetics:</b>	17 B1

**Gene Summary:**

Mediates interactions of advanced glycosylation end products (AGE). These are nonenzymatically glycosylated proteins which accumulate in vascular tissue in aging and at an accelerated rate in diabetes. Acts as a mediator of both acute and chronic vascular inflammation in conditions such as atherosclerosis and in particular as a complication of diabetes. AGE/RAGE signaling plays an important role in regulating the production/expression of TNF-alpha, oxidative stress, and endothelial dysfunction in type 2 diabetes. Interaction with S100A12 on endothelium, mononuclear phagocytes, and lymphocytes triggers cellular activation, with generation of key proinflammatory mediators. Interaction with S100B after myocardial infarction may play a role in myocyte apoptosis by activating ERK1/2 and p53/TP53 signaling. Can also bind oligonucleotides. Receptor for amyloid beta peptide. Contributes to the translocation of amyloid-beta peptide (ABPP) across the cell membrane from the extracellular to the intracellular space in cortical neurons. ABPP-initiated RAGE signaling, especially stimulation of p38 mitogen-activated protein kinase (MAPK), has the capacity to drive a transport system delivering ABPP as a complex with RAGE to the intraneuronal space. RAGE-dependent signaling in microglia contributes to neuroinflammation, amyloid accumulation, and impaired learning/memory in a mouse model of Alzheimer disease. [UniProtKB/Swiss-Prot Function]

Transcript Variant: This variant (4) contains multiple differences in the coding region, compared to variant 1, and initiates translation at an alternate start codon. The encoded isoform (d) is shorter and has a distinct N-terminus, compared to isoform a. Sequence Note: The RefSeq transcript and protein were derived from genomic sequence to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on alignments.