

## Product datasheet for **MC227001**

### Ankrd60 (NM\_027303) Mouse Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Ankrd60 (NM_027303) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Ankrd60
Synonyms:	1700019A24Rik; 1700030G11Rik; RP23-479E23.1
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
Fully Sequenced ORF:	>MC227001 representing NM_027303 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**GCGATCGCC**

ATGACGCGGGGATGACGCGGGCCATGCGGAGGGCGGCAGACGTGGCGGGGGCGGGCGCGCTCCGCCGA  
AACGGCAGGGGGACGACCCCAAGGCGGGCGCGCCACACTCAGCCCTCCAGCAGGGCCCCGCGCACAGGT  
TAGGGGGCACACAGGGGACCCGCTTGCCCTCCACCGTGTGCCAGACTTGACTCCTGACGTCTATGTCCTG  
CGGGTGAAGCTGGAGGAGACTGGGGAGCTGTCCGCGTGGCCAACTGCCGCAACGACATGACAGTGCGGG  
AGCTCAAAGAAGAGCTGGACCTGATGGTCGGCATCCCCCTCAACCTGCAGCGGCTGCATTTCTGGACCG  
AGGAATTCGTATGGATGATGCTACCCTGAAGTTCTACGATGTCATCCCCGGAGCACTTATTTATTATGT  
ATCTGGCATTACGATGGTTGGACGGAAGTGGTTTTGGCCGCCGTGGAAGGGGACCCAGCAAGATGTCAT  
GTCTTGGCATCAGCGAGGATACCTTCTACCGAACAGGCAACTCACAGCATTTTAAGGGGGAACAGTGGAG  
GCAGTGGACAGCCAGAGGGCATTGTGGCCTGTACATCTCTGCTCACAGAGGCCACGTGGATGCCGTG  
CAGTATCTTAGGACTTGGAGCAAACCTGCCTGGCCAAATCACCCCTGGGCAGGACGCCCTGCACGTGG  
TCGCAGCCATGGGCAGGACAGACTGCGTCCAGCCCTTACTGGAACATGGTGCCTCCATCCATGACAGAGA  
CTTGAAGAGCGACACGCCATCACCATTTGCGCGGCGCCTGAAGCGAGGACATTTAGAGCACAAAGATGTTT  
CTTCTCTACTGGATGGTCAAGTCGGGGTCCCGGGACCCGACTTGGTAGGGAGGAAAGCTTTGGAGAAAT  
CCAGCCCTGGGATGGCGACCAGGAGGAGACAATATTAATCTCCTGA

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	Sgfl-MluI
ACCN:	NM_027303
Insert Size:	957 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_027303.1</a></u> , <u><a href="#">NP_081579.1</a></u>
<b>RefSeq Size:</b>	1065 bp
<b>RefSeq ORF:</b>	957 bp
<b>Locus ID:</b>	70065
<b>Cytogenetics:</b>	2 H3