

## Product datasheet for **MC227017**

### Gm6408 (NM\_001243104) Mouse Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Gm6408 (NM\_001243104) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Gm6408  
**Synonyms:** 623198  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >MC227017 representing NM\_001243104  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**GCGATCGCC**

ATGTTCTCCTGTTTCCAGGGTTCAGAGGCTCCGGCCACAAAAAGCCAAGAGCGGTTTCCTTGTTTCGCT  
TTTGGAGACGTTTGATTCGCCCTCTGACACTTCAGGCATGCTCCACAGCGAGCCCAAGGTCTGTTT  
TCAAAATGAACAGGAGCCAGACTGTATGCCAGCGCCCGGATTCAACTATAACAGCCAGAATATGTG  
GTACAAATGATCCACTACATTCCAGCTGCTGTTTCAAGCATAATGACCATGTCTGTATCCGCATTTTCTTG  
CCATGTACCCATCGTACGCTAGCACTTGAAGGTGCTGGACCTGCTGATGACAACGTATGCATTCTTCT  
GCCTGACCGTATAAAGGACCAAGAACTAAGAGAGCCATATTCAGGTTTCTATTCCACTGGTTTAAAGAAA  
TTCCCCAAGGATTTTATGAGTCCCCAGACCTGGCTGTTGTGAGGCAGTTCATTGACTATGTGAGGCATA  
ATGTGCCTTCCGACAGCAAGATACCCAAGCCAGGGAGCTGCTCTCTGTGCTGGAGGAGCAGGAGGCCAT  
TGGGCTAAACCTGAGGAAGATTTTCCACAGCTCCAGAGCCCTCAGAGCAGGATGCCTCTGGGAGGCAG  
GAGACTTGGCTCCGAGGCCAGCTTCTGAGACAGAGCCACAGGGAGACAAGCAACCCCGAGAGGATCCTG  
AGCTTGTGAAAGGAGCAGTTCTGGATCCCTCTGAACCAAAAGCCACGATAGAGCTGCCTCCGTCTTTACA  
TCAAGCTGTGCCACATCAGATGGTACACTCTACCTGTAGATGTTACTGCAGATGAGGCGACAGGTGTG  
GCTGCAGATGAGGCTGCAGATGTTGCTGCAGATGTTTACCTGCCAGGCAGTTATTTGTGCTTACACTG  
TGCAGTTAGTACACCAGATTTTGTATCCCTCTGCCTGAGGTTGACAT**ATAG**

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_001243104  
**Insert Size:** 963 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001243104.1</a></u> , <u><a href="#">NP_001230033.1</a></u>
<b>RefSeq Size:</b>	1191 bp
<b>RefSeq ORF:</b>	963 bp
<b>Locus ID:</b>	623198
<b>Cytogenetics:</b>	5