

## Product datasheet for **MC227059**

### Cdk19 (NM\_001291817) Mouse Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Cdk19 (NM\_001291817) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Cdk19  
**Synonyms:** 2700084L06Rik; AI845279; AW228747; Cdc2l6; CDK11; mKIAA1028  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >MC227059 representing NM\_001291817  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**GCGATCGCC**

ATGGGTTTTGCCAGGTTATTCAATTCTCCCCTAAAGCCACTCGCAGATTTGGATCCAGTGGTTGTGACAT  
TTTGGTATCGGGCTCCGGAACCTTTACTTGGTCCAGGCATTACACAAAGGCCATTGACATCTGGGCAAT  
AGGCTGCATATTTGCTGAACTCTTGACTTCAGAACCCATTTTTACTGTGCTCAGGAGGATATAAAAACA  
AGCAATCCTTTTCATCATGATCAGTTAGATCGAATATTTAGTGTCATGGGGTTTCTGCAGATAAAGACT  
GGGAAGATATTAGAAAAATGCCAGAGTACCCAACACTTCAGAAAGACTTTGAAAGAAACAGTACGCCAA  
CAGCAGCCTCATAAAAATACATGGAGAAGCACAAGGTCAAGCTGACAGCAAAGTGTCTCTGCTTCAG  
AACTCCTCACCATGGATCCAACCAAGAGAATCACCTCAGAGCAGGCTCTGCAGGACCCTTACTTCCAGG  
AAGACCCCTGCCAACATTAGATGTGTTTCGCTGGCTGCCAGATTCGTAACCCAAAGAGAATTCCTTAA  
TGAAGATGAACCAAGAGAAAGGTGACAAGAACCAGCCACAGCAGCAGAACCCACATCAGCAACCTGCA  
GCCCCGCACAGCAGACAGCAGCCCCCACAGGCTCCTCCACCACAGCAGCAGTGCACCAACCAATG  
GAACTGCTGGGGAGCCACTGCGGGAGGGGGGGCGCTGGAGCGGGCTGCAGCACAGCAGGACCCCGG  
GCTGAATCAGGTGCCTCAAACAAGAAGCCCGGATTGGGCTTCCGGTGCAAACTCAGGCGGGCCGTT  
ATGCCGTCGGATTATCAGCATTCCAGTTCTCGCCTGAATTACCAAAGCAGCGTACAGGGTCTCTCAGT  
CCCAGAGTACGCTAGGCTACTCTTCATCTCAACAGAGTACCCAGTACCACTCATCTCACCAGACCCACCG  
GTACTGA

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_001291817



[View online »](#)

<b>Insert Size:</b>	987 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001291817.1</a></u> , <u><a href="#">NP_001278746.1</a></u>
<b>RefSeq Size:</b>	5770 bp
<b>RefSeq ORF:</b>	987 bp
<b>Locus ID:</b>	78334
<b>Cytogenetics:</b>	10 B1