

Product datasheet for **MC227073**

Tmem255a (NM_001289727) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Tmem255a (NM_001289727) Mouse Untagged Clone
Tag: Tag Free
Symbol: Tmem255a
Synonyms: 4933417N17; 6430550H21Rik; Fam70a
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC227073 representing NM_001289727
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGCATCAGTCCCTGACTCAGCAGCGGTCCAGCGACATGTCCTTCCCGATTCCATGGGAGCCTTCAATA
GGAGGAAGCGGAAGTCCATCTACGTCACCGTTACTTTGCTTATTGTGTCATGTTAATTCTCACGGTGGG
CCTTGCTGCGACTACCAGGACCCAGAATGTGACCGTGGGAGGTTATTACCCAGGAGTTATTCTTGCTTT
GGATCATTCCCTGGGATCATCGGGTCAAACCTGATTGAGAACAAGCGGCAGATGCTGGTGGCTTCTATCG
TGTTTATCAGCTTTGGTGTGATTGCAGCATTTTGTGTGCTATAGTTGATGGGGTCTTGGCCGACAGACA
CATTGATTTGAAACCACTCTATGCTAACCAGGTCGCACTATGTTCCCAAGACATCAGAGGGAAGCTGAG
GAGGTCATAAATTCCCTATCCAAAATTACTCCTTCTACGAGGGCTTTGAGGAACCTTACCAGGCAGTTA
AAGAGGTGAAGTGCCTCAACTTAGCCGTGGACTTTGCACACCTCGCATCCGGGGCAACACTGCTTCTG
CTGTGACCTCTACAAGTGGCAACAGGGTGGAAATCACCGGCGGCTACTATGAGTACATTGATGTCAGC
AGCTGCCAGGACATCATCCACCTGTACCACCTACTCTGGTCTGCCACCCTCAACATCGTTGGCCTGT
TCCTGGGTATCATTACTGCTGCGGTCTTGGGGATTTAAGGATATGAACCTACTTCCGGCACTGAA
CTGCTCTGTTGAAAATGCCATCCAAGTCTTCTTACTATGCGCGTCCCAAGTGGCATCGTACAATACC
TACTACCATAGCCCTCCTCATCTGCCACCGTATTCTGCTTATGACTTTACGGTATTCTCGTGGTGTGTTG
AAAAGTACAGCAGGCACTCTAAAACAAGTGGCAATGCCAAGTACTTTAACTGCAATGGTTCAAGAAAC
AATGGGACGTGCT**TAA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001289727



[View online »](#)

| | |
|-------------------------------|---|
| Insert Size: | 996 bp |
| OTI Disclaimer: | Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP). |
| OTI Annotation: | Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag. |
| Components: | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water). |
| Reconstitution Method: | <ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C. |
| RefSeq: | <u>NM_001289727.1</u> , <u>NP_001276656.1</u> |
| RefSeq Size: | 1797 bp |
| RefSeq ORF: | 996 bp |
| Locus ID: | 245386 |
| UniProt ID: | <u>Q8BHW5</u> |
| Cytogenetics: | X A3.3 |