

Product datasheet for **MC227081**

Fam131b (NM_001286584) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Tag Free
Symbol:	Fam131b
Synonyms:	6330503C03Rik; 6530406I18Rik; AV277466
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>MC227081 representing NM_001286584 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGGACAGCACCAGCTCGCTGCATGGCAGCAGTCTCCATCGGCCTTCGACAGAGCAAACACGAACAGATT
TCTCCTGGGACGGCATCAATCTTTCCATGGAGGACACCACTTCCATCCTCCGAAGCTTAAGCGAAACT
TAATGCCTATGGCATTGGGGCCCTGGCCAAGTCATCATTCTCAGGGATCTCGAGGAGCATGAAGGACCAT
GTGACAAAGCCTACCGCATGGGCCAAGGCCGGTGGCCACATGATTGAGTGGCAGGGCTGGGGGAAGG
CTCCAAACAATTCAGCCACAACATAGCCATGAGGCTGTGCGCAGGGATACAGACGCCTACTCTGACCTCAG
TGATGGCGAGAAGGAGGCACGTTTCTAGCAGGTGTCATGGAACAGTTTGTATTTCTGAGGCTACACTT
ATGGCTTGGTCTTCCATGGATGGTGAAGACATGAGTGTCAATTCTACCCAGGAGCCACTGGACTGCAACT
ACAGTGACAACCTACCAGGAGCTGATGGAAAGTCAGGATGCACTAGCTCAAGCCCAATGGATGGATGGCC
TCACTCCTACGTGTACAGGGGATGTACTGTCTGGGGTCATCTGACGCCTGGGAAGCAAGTGATCAGTCC
CTCATTGCATCTCCAGCTACCGGATCCTATCTTGGCCCTGCATTTGATGACTCACAGCCGAGCCTGCATG
ACATGGGGCCTTCCAACCTGCTTCAGGATACTCTGCTCAGGAGCCTCCACCTTTGTTGGGGGTAGACAC
TGACTGGGCATCAGAGGTTGGTGGGGTGGAACTGGCCAGGGGACCTGTAGAGGAAGAGAAGAGACCCTG
GCCCTGAGGAGGAAGAAGATGCAGGATGCCGGGATCTGGAGTCACTTTCCCCACGAGAGGCCCAAGAGA
TGTCCACTGCTCTCAGCCGGAAGGTGTCTGATGTCACATCCTCAGGTGTGCACTCTTTGATGAGGAGGA
AGGGGATGCCAACAAC**TAG**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI



ACCN:	NM_001286584
Insert Size:	999 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_001286584.1</u> , <u>NP_001273513.1</u>
RefSeq Size:	4048 bp
RefSeq ORF:	999 bp
Locus ID:	76156
UniProt ID:	<u>Q3TY60</u>
Cytogenetics:	6