

Product datasheet for **MC227087**

Gpr146 (NM_030258) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Gpr146 (NM_030258) Mouse Untagged Clone
Tag: Tag Free
Symbol: Gpr146
Synonyms: BC003323; PGR8
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Fully Sequenced ORF: >MC227087 representing NM_030258
Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTGGAGCTGTGGCCCACTCAACAGCACAGCGTGGGCTGAGGAGCCGCTGTGCCGAACCTGCGCCTGG
GGCTGTGGGTCTCTCGTGTCTACCTGGGGCAGGGTCCCTGTGAGCTTAGGCTATAATGCTTTCT
GGTGTGGCAACCTGGCCAGCAAGAACCATTGACCATGCCCGACGTGACTTCGTGAACATGGCCGTG
GCGGGGCTGGTGTCTACGGCACTGGCACCTGCGTACCTGTTGGGCCCTGCCCACTCCAGGTGGGCCCTGT
GGAGCCTCAGCAGTGAGGCCCATGTGACACTGCTCATCCTGTTCAACGTGGCTTCCTGGTGACCATGTA
CTCCACTGCACTGCTGAGCCTTGACTACTACATCGAGCGTGCCCTGCCACGCACCTACATGGCCAGTGTG
TACAACACCCGGCACGTGTGTGGCTTCGTCTGGGAGGGGCGGTGCTCACCAGCTTCTCCTCCCTGCTCT
TCTACATCTGCAGTCACGTGTCTTCTAGAATCGCTGAGTGTGCCCGGATGCAGAACACGGAGGCAGCCGA
TGCTATCCTTGTGCTCATCGGCTACGTGGTGCCAGGTCTGGCTGTGTTGTATGCCCTGGCACTCATCTCG
AGAATCGGGAAGGAAGACACACCCTGGACCAGGACACAGCAGGCTGGACCCCTCGGTGCACAGGCTGC
TGGTGGCCACCGTGTGCACTCAGTTTGGCCTCTGGACACCTTACTACTTGAGCCTGGGGCACACAGTGCT
GACGTACGCGGGGAGGACCGTGGAGGGGCATTATCTGGGCATCTACAGTTGCTAAGGACCTGGCCAAG
TTCCTGGCCTTCTCAAGCAGTTCTGTGACCACTGCTCTACCGTTACATCAACAAAGCCTTCCCCGGCA
AGCTCCGGCGGCTGATGAAGAAGATGCACTGCGGGCGCCGCACTGCTCCCCGACCCCTCGGGGATACA
GCAGGTGATGGCACAGGCG**TAG**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI



[View online »](#)

ACCN:	NM_030258
Insert Size:	1002 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_030258.4</u> , <u>NP_084534.2</u>
RefSeq Size:	4065 bp
RefSeq ORF:	1002 bp
Locus ID:	80290
UniProt ID:	<u>Q99LE2</u>
Cytogenetics:	5 G2
Gene Summary:	Orphan receptor.[UniProtKB/Swiss-Prot Function]