

Product datasheet for **MC227088**

Cryz11 (NM_026994) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Cryz11 (NM_026994) Mouse Untagged Clone
Tag: Tag Free
Symbol: Cryz11
Synonyms: 2210407J23Rik; 2410006O11Rik; QOH-1
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Fully Sequenced ORF: >MC227088 representing NM_026994
Red=Cloning site **Blue**=ORF **Orange**=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGAAAGGCTTGTATTTTCAACAGAGTTCTACAAATGAAGAAGTAAACATTTGTATTTCAAGAAAAGGAAA
ATGTCCCTGTCACAGAGGATAACTTTGTGAGAGTTCAAGTTAAAGCTTGTGCCCTGAGCCACATAAACAC
AAAGCTTCTAGCAGAAAATGAAGATGGAAAAGGATTTCTCCCTGTTGGAAGAGAAGTTTCTGGAATTGTA
TTAGAGGGAATTCTGCCTTTGGATTCCGAAGACCCCGGGCTCTGTGAGGTTATTCGAGTACATGAGCATT
ACCTGGTTCATAAGCCAGAAAAGGTGTCATGGACGGAAGCTGCTGGCGTCATCCGGGACGGAGTCCGAGC
CTGCACAGCTCTGTATTATCTTTCTCAACTCTCTCCTGAAAGTCTGTGCTGATCATGGACGGGGCAAGT
GCATTTGGCACAATAGCCATTAGTTAGCCACCACAGAGGAGCCAAGGTTATCTCAACAGCACACAGCC
TGGAGGACAAGCAACATCTTGAAGGCTTAGGCCCTTATAGCCCGTGTGATTGATGTGTCTAATGGGAA
AGTCCATGTTGCTGAAAGCTGTTTGGAGGAAACAGGTGGCCTGGGAGTAGATATTGTGATAGATGCAGGA
GTGAGATTATATAGTAAAGATGATGAACCAGCTGTAAAGTTACATCTGCCACATAAGCATGACATCATCA
CACTTCTTGGTGTCCGAGGCCATTGGGTAACAACAGAAGAAAACCTTCAGTTGGATCCTCCAGATAGCCA
CTGCCTTTTCTCAAGGGAGCAACGGTGGCTTTCTTAATGATGAAGTTTGAATTTGTCAAATGCACAA
CAGGAAAAATATCTTTGTATCTTAAAAGATGTGATGGAGAAGTTATCAGCTGGTGTTCAGGCCTCTGC
TCGATGAGCCCATCCCACTGTATGAGGCAAGGTTTCCATGGAAGTTGTTGAGAAAAATCAAGAAAAGAAA
AAAGCAAGTCGTTCAATTT**TAA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI



ACCN:	NM_026994
Insert Size:	1002 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_026994.1, NP_081270.1</u>
RefSeq Size:	1747 bp
RefSeq ORF:	1002 bp
Locus ID:	66609
Cytogenetics:	16 C4