

Product datasheet for **MC227090**

Mavs (NM_001206383) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Mavs (NM_001206383) Mouse Untagged Clone
Tag: Tag Free
Symbol: Mavs
Synonyms: cardif; D430028G21Rik; IPS-1; Visa
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC227090 representing NM_001206383
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGTCTGGTGGCTCTTTGATACCTCTCCTAACCAGCAGGCTCTCAGCCCTCAGCCCTCCAGAGAGCATC
 AAGAGCAAGAACCAGAAGTGGTGGCGCCACGCAGCAAATGTTGCCTCTGTTCCCATAGCAACCTATGG
 ACCTGTGTCTCCAACCGTTTCCTCCAGCCCTTCCACGTAAGTGCCTGAGGACAAACCTTTGTCTGGG
 GTCACAGTATCAGCCCTATCTGCTGATACCTCTTTGTCCTCCTCGTCCACTGGATCAGCTTTGCAAAGG
 GAGCTGGTGACCAGGCCAAAGCTGCCACCTGTTTCAGTACTACACTCACCAATTCTGTGACTACCGCTC
 AGTGCCTTCTCCAGATTGGTCCAGTAAAAACCATGTCTTCCAAGTTGCCCTCAGTTCAAAGTCCACT
 GCTGCGATGACGTCTACTGTGCTACCAATACAGCGCCATCAAAATTACCCAGCAACTCAGTGTATGCGG
 GCACAGTGCCATCCAGAGTGCCTGCTAGTGTGGCCAAAGCACCTGCCAACACAATACCACCTGAGAGGAA
 CAGCAAGCAAGCCAAGGAGACCCCGAGGGTCCAGCAACCAAGTCACTGGAGGCAACCAGACTGGA
 CCAAATAGCAGTATCAGGAGCTTGCACTCTGGACCAGAGATGAGCAAGCCAGGTGTGCTGATCCAGT
 TGGACGAGCCATTCTCAGCCTGCTGTGGACCTTGCCATTAGCCCTAGCAGCTCCTTGGTCTCAGAAC
 CAACCATGGTCCAGAGGAGAATGAGTATTCGTCCTTTAGAATCCAGGTAGACGAAAGCCCAAGTGCCTGAT
 CTATTAGGAAGCCCTGAGCCACTAGCCACCCAGCAGCCCAAGAAGAGGAAGAACATTGTGCCAGTTCAA
 TGCCCTGGGCTAAGTGGCTTGGGGCCACCAGTGCCTTGGCTGTATTCTGGCAGTGTGCTGTACCG
 TAGTAGGCGCCTGGCCAG**TGA**

ACGGTACGGCGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001206383



[View online »](#)

Insert Size:	1002 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001206383.1</u> , <u>NP_001193312.1</u>
RefSeq Size:	2696 bp
RefSeq ORF:	1002 bp
Locus ID:	228607
Cytogenetics:	2 F1
Gene Summary:	<p>Required for innate immune defense against viruses (PubMed:24037184). Acts downstream of DHX33, DDX58/RIG-I and IFIH1/MDA5, which detect intracellular dsRNA produced during viral replication, to coordinate pathways leading to the activation of NF-kappa-B, IRF3 and IRF7, and to the subsequent induction of antiviral cytokines such as IFN-beta and RANTES (CCL5) (PubMed:24037184). Peroxisomal and mitochondrial MAVS act sequentially to create an antiviral cellular state (By similarity). Upon viral infection, peroxisomal MAVS induces the rapid interferon-independent expression of defense factors that provide short-term protection, whereas mitochondrial MAVS activates an interferon-dependent signaling pathway with delayed kinetics, which amplifies and stabilizes the antiviral response (By similarity). May activate the same pathways following detection of extracellular dsRNA by TLR3 (By similarity). May protect cells from apoptosis (By similarity).[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (4) differs in the 5' UTR, lacks a portion of the 5' coding region and initiates translation at a downstream, in-frame start codon, compared to variant 1. Variants 3 and 4 encode the same isoform (2), which has a shorter N-terminus, compared to isoform 1. Sequence Note: The RefSeq transcript and protein were derived from genomic sequence to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on alignments.</p>