

Product datasheet for **MC227105**

Gpr160 (NM_001286995) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Gpr160 (NM_001286995) Mouse Untagged Clone
Tag: Tag Free
Symbol: Gpr160
Synonyms: 1700025D19Rik; GPCR150
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC227105 representing NM_001286995
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**GCGATCGCC**

ATGACAGCTCTTTCTTCAAAGAAGTCTCTCTTCAATACCAGCTCCATCAGTCTCCGAGCTCCTGGAGG
 CCTCCTGTCTGTTCTTGATCATCCTGGGAAAGTGTATTAACATCCTCCTCTAAGAGTGAGAAG
 GGGAGACGCCCGCTGGACTTTGATGGAGTACTTCTGCTTCTCACTCGCACTTGTGATCTCCTGCTGTTG
 GTAACATTTCCATCCTGACCTACTTCAGGGACTTTGTAGTTCTAGGTATTAGGTTTACCAGATACCACA
 TCTGTCTGCTCACACAAATTATTTCTTCACTTATGGCTTCTTGCATTACCCAGTCTGCTCACTGGCTTG
 TATAGACTACTGGTGAATCTGTCCAGAGCCAGCAAGCAGTCACTAGGTGGCAGAAATTAATTTATTTTC
 TTGACAGTCATTTAACTTGGATCTCAGTCCTTGCTTACGTTTTGGTTGACCCGGCCATCTCAGTCAGCC
 TAAAGGCACACAGGGGCTATGTGTACCAGTGCCCGCCTATGTCAGCACTCAGAGTCACTGGCTCTCACT
 GTCTATGCTGATGGTTTTATTTGTGGCTTTTCTGATTTTCATGGCAGGAAGTTGTTGCCTTGCTTCAGGT
 ATGAGGATAGCATCCTATAAGAGCAAGGCCGCCCTACTTTCCCTTCCCACTCCACTGCCGTTACGCC
 TGAGCTGCAGAGAAGCACTTTGCCAGGCTCATCGTGTGCTTTCTTGGGACCTGGTTCCCTTTGTAGC
 ACTTCAGTCTCATCCTCTCCCTCAGAGTTCAGATCCCGCGGTACATAGAAATGAACGTCCCTGGCTG
 TACTTTGTCAACAGTTTCTCATCGCTGCAGTTTATTGGTTTAACTGCCACAAGCTTGATTTAAGAGACA
 GCTCGTTACCTGTGGATCCTTTCACTCAACTGAAATGCTGCTTTGTTCCAGTCCATAGACTTAAGCAAGT
 GGAGAGGCCCATGTCCATAGTCATCTGT**TAG**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001286995



[View online »](#)

Insert Size:	1011 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001286995.1</u> , <u>NP_001273924.1</u>
RefSeq Size:	1920 bp
RefSeq ORF:	1011 bp
Locus ID:	71862
UniProt ID:	<u>Q3U3F9</u>
Cytogenetics:	3 A3
Gene Summary:	Orphan receptor.[UniProtKB/Swiss-Prot Function]