

## Product datasheet for **MC227117**

### **P2ry14 (NM\_001287119) Mouse Untagged Clone**

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** P2ry14 (NM\_001287119) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** P2ry14  
**Synonyms:** A330108O13Rik; Gpr105; P2Y14  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >MC227117 representing NM\_001287119  
**Red**=Cloning site **Blue**=ORF **Orange**=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**GCGATCGCC**

ATGAACAACCTCCACCACCAGACCCCTCAAACCAGCCCTGCTCCTGGAACACCCTGATCACAAGCAGA  
TCATCCCGTGTGTACGGTATGGTCTTCATCACGGGGCTCCTCCTCAATGGGATATCAGGATGGATATT  
CTTTTATGTGCCAGCTCCAAGAGTTTTATCATCTATCTCAAGAACATAGTGGTGGCTGACTTTCTCATG  
GGCCTGACTTTCCCTTTCAAAGTCCTTGGTGACTCAGGCCTCGGCCCTGGCAGGTGAATGTGTTTGTGT  
GCAGGGTCTCTGCCGTATCTTCTATGTTAATATGTACGTCAGCATCGTGTCTTTGGGCTCATCAGCTT  
TGACAGGTAATAAAAATTGTGAAGCCCTTCTGACGCTATTGTGACGTCGGTGAATATAGCAAGCTG  
CTTTCTGTGCTCGTGTGGATGCTCATGCTTCTCCTTGCTGTCCAAACATCATCCTGACAAACCAGGGTG  
TCAAGGAGGTCACGAAGATACAGTGCATGGAGCTCAAAAACGAGCTGGGGCGGAAGTGGCACAAGGCGTC  
TAACTATATCTTCGTGAGTATCTTCTGGGTCGTGTTTCTTCTGCTAATCGTCTTCTACACGGCCATCAGC  
AGGAAGATCTTCAAGTCTCACCTCAAGTCCAGGAAGAATCCACCTCCGTCAAGAGGAAGTCCAGCCGCA  
ATATCTTCAGCATCGTGCTCGTTTTTGTGCTCTGCTTTGTGCCTTACCACATTGCCAGAATCCCTTACAC  
AAAGAGTCAGACGGAAGTCACTACAGCTGCCGACGAAGGAGACCCTGCTCTATGCGAAAGAATTCATC  
CTGCTACTCTCGGCTGCCAATGTGTGTCTGGACCCATTATTTATTTCTTATGCCAGCCATTTAGAG  
AAGTCTTAAATAAGAAGTTACACATGCTCAAGTCCAGAATGACCTAGAGGTTTCCAAAACAAAAG  
GGAAAATGCGATTATGAAAGCACAGATACTTTGTAA

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_001287119



[View online »](#)

<b>Insert Size:</b>	1017 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001287119.1</a></u> , <u><a href="#">NP_001274048.1</a></u>
<b>RefSeq Size:</b>	1617 bp
<b>RefSeq ORF:</b>	1017 bp
<b>Locus ID:</b>	140795
<b>UniProt ID:</b>	<u><a href="#">Q9ESG6</a></u>
<b>Cytogenetics:</b>	3 28.96 cM
<b>Gene Summary:</b>	Receptor for UDP-glucose coupled to G-proteins.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (3) differs in the 5' UTR compared to variant 1. Variants 1 through 7 encode the same protein (isoform a). Sequence Note: The RefSeq transcript and protein were derived from genomic sequence to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on alignments.