

Product datasheet for **MC227118**

P2ry14 (NM_001287120) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: P2ry14 (NM_001287120) Mouse Untagged Clone
Tag: Tag Free
Symbol: P2ry14
Synonyms: A330108O13Rik; Gpr105; P2Y14
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Fully Sequenced ORF: >MC227118 representing NM_001287120
Red=Cloning site **Blue**=ORF **Orange**=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGAACAACCTCCACCACCACAGACCCTCCAAACCAGCCCTGCTCCTGGAACACCCTGATCACAAAGCAGA
TCATTCCTCGTGTGTACGGTATGGTCTTCATCACGGGGCTCCTCCTCAATGGGATATCAGGATGGATATT
CTTTTATGTGCCAGCTCCAAGAGTTTTATCATCTATCTCAAGAACATAGTGGTGGCTGACTTTCTCATG
GGCCTGACTTTCCCTTCAAAGTCCTTGGTGACTCAGGCCTCGGCCCTGGCAGGTGAATGTGTTGTGT
GCAGGGTCTCTGCCGTACTTCTATGTTAATATGTACGTCAGCATCGTGTCTTTGGGCTCATCAGCTT
TGACAGGTAATAAAAATTGTGAAGCCCTTCTGACGCTATTGTGCAGTCGGTGAATAAGCAAGCTG
CTTTCTGTGCTCGTGTGGATGCTCATGCTTCTCCTTGCTGTCCAAACATCATCCTGACAAACCAGGGTG
TCAAGGAGGTCACGAAGATACAGTGCATGGAGCTCAAAAACGAGCTGGGGCGGAAGTGGCACAAGGCGTC
TAACTATATCTTCGTGAGTATCTTCTGGGTCGTGTTCTTCTGCTAATCGTCTTACACGGCCATCACG
AGGAAGATCTTCAAGTCTCACCTCAAGTCCAGGAAGAATTCCACCTCCGTCAGAGGAAGTCCAGCCGCA
ATATCTTCAGCATCGTGCCTGTTTTGTGCTGCTTTGTGCCTTACCACATTGCCAGAATCCCTACAC
AAAGAGTCAGACGAAGGTCACCTACAGCTGCCGGACGAAGGAGACCCTGCTCTATGCCAAGAATTCACT
CTGCTACTCTCGGCTGCCAATGTGTGCTGGACCCATTATTTATTTCTTCTTATGCCAGCCATTTAGAG
AAGTCTTAAATAAGAAGTTACACATGCTCACTCAAAGTCCAGAATGACCTAGAGGTTTCCAAAACAAAAG
GGAAAATGCGATTTCATGAAAGCACAGATACTTTGTAA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI



[View online »](#)

ACCN:	NM_001287120
Insert Size:	1017 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001287120.1, NP_001274049.1</u>
RefSeq Size:	1742 bp
RefSeq ORF:	1017 bp
Locus ID:	140795
UniProt ID:	<u>Q9ESG6</u>
Cytogenetics:	3 28.96 cM
Gene Summary:	Receptor for UDP-glucose coupled to G-proteins.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (4) differs in the 5' UTR compared to variant 1. Variants 1 through 7 encode the same protein (isoform a).