

Product datasheet for MC227151

Cd2bp2 (NM_001285907) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Cd2bp2 (NM_001285907) Mouse Untagged Clone
Tag: Tag Free
Symbol: Cd2bp2
Synonyms: 1500011B02Rik; 2410024K20Rik; AW048183
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC227151 representing NM_001285907
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGC**C

ATGCCAAAGAGGAAAGTGACTTTTCAGGGTGTGGAGATGAGGATGGTGAAGATGAAATCAGTGTGCCA
 AAAAGAAGCTGGTGGACCCTGTGGCTGCAGCAGGTGGTCTGGGAGCCGCTTCAAAGGCAAGCACTCTTT
 AGACAGTGATGAAGAGGATGATGATGAGGAAGGGTCCAGCAAATATGATATCCTGGCCTCAGAGGATGTA
 GAAGGTACAGGAAGCAGCTACTCTCCCTAGTGAGGGCGGTGTGAGGATCACACCCTTCAACCTGCAGGAAG
 AGATGGAGGAAGGCCACTTTGATGCTGATGGCAACTACTTCTGAACCAGGATGCCAGATCCGAGACAG
 CTGGCTGGACAACATTGACTGGGTGAGGATCAAGGAACGGCCCCCTGATAAGCACCAGGTTTCAGACTCG
 GAGGAGGAAGACAGCCTGGGTGAGACTCCAATGAGCGCCCAAGCACTCCTGGAGGGACTTTTGGAGTTGT
 TGTGCCAAGAGAGACAGTGGCTGGGGCCCTTAGGCGGCTGGGGCCCGAGGAGGAGGCAAAGGGGAGCAA
 CAGCAAGGGGACTGGACGGCCCAATCCCCCAGCGGCTAGACCGGCTGTCTGGTTAGCTGACCAATG
 GTGGCACGAGGCAACCTTGGTGTGTACCAGGAAACAAGGAACGGTTGGCTATGCGACTGAAGGGTTGG
 GGTGTCCGGCCAGGGATCCCATGACCCACACCTCCACCTCCCTGGACATGTTTGTGAGGAAGTAGC
 AGAGGGTGAACCTGAAACCCCAACCCCTACCCAAAGAGAGGAGGCAGAGTCAGCAGGAGATGGCCTGATG
 GATGTGATGTGGGAGTATAAATGGGAGAACACGGGGACGCCGACTGTATGGACCCTTACCAGTGCC
 AGATGCAGACGTGGGTGAGTGAAGGTTACTTCCGGATGGTGTATTGTGCGAAGCTGGACCCTCCGGG
 CGGCCAGTTCTACAACCTCAAGCGCATTGATTTTGGAGCTGTACACG**TGA**

ACGGTACGGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001285907



[View online »](#)

Insert Size:	1029 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001285907.1</u> , <u>NP_001272836.1</u>
RefSeq Size:	3509 bp
RefSeq ORF:	1029 bp
Locus ID:	70233
UniProt ID:	<u>Q9CWK3</u>
Cytogenetics:	7 F3
Gene Summary:	<p>Involved in pre-mRNA splicing as component of the U5 snRNP complex that is involved in spliceosome assembly.[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (1) represents the longest transcript. Variants 1-5 encode the same protein. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>