

Product datasheet for **MC227159**

B4galt4 (NM_001285793) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: B4galt4 (NM_001285793) Mouse Untagged Clone
Tag: Tag Free
Symbol: B4galt4
Synonyms: 9130402O08Rik; b4Gal-T4; B4galt-IV; beta4Gal-T4; beta4GalT-IV
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC227159 representing NM_001285793
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGGCTGCAACCCACCTTATCACCTCTCCTACAGATTACGATTGTTGCTGCTTTACCTGTGCCTGA
 CGGTGGTGGGTGGCCACCAGCAACTACTTTGGGTGCTATCAAGTGATCCCCAAGGCAAAGGACTT
 CATGGCTAGTTTCCACAAGGTCATACATTTGGGAATGAAGAACTCTGGCCATGATGGGGCCACGAAA
 AAACCAGAGCTTGCTAACTGCCCTTCGGTGTCTCCAACTCAGAGGCCAGAGCAAGCTCGTTTTAAGC
 CGGACCTCACGTTGGAGGAAATAGAGGCTGAAAACCCCAAAGTGTCCAGAGGCCGGTATCACCTGAGGA
 ATGTAAGGCTCTGCAGCGGGTGGCCATCCTCATTCCACACAGGAACAGAGAGAAGCACTTGATCTACCTG
 CTGGAACACCTGCATCCCTTCCTGCAGAGGCAACAGCTGGACTACGGCATCTATATCATCCACCAGACGG
 GAAGTAAAAAGTTTAAACCGAGCCAAGCTCCTGAACGTGGGCTATCTGGAAGCTCTCAAGGAGGAACTG
 GGACTGCTTCGTATCCACGACGTGGACCTGGTGCCTGAGAATGACTTCAACCTCTACACCTGCGGTGAT
 CAGCCCAAGCACTTGGTGGTGGCCGGAACAGCACGGGCTACAGGTTGCGTTACAGTAAATATTTGGGG
 GTGTCCTGCCCTCAGCAGGGAACAGTTTCTCAAGGTGAATGGATTCTCTAACAACTACTGGGGTGGGG
 AGGAGAAGACGATGACCTCAGACTCAGGGTTGAGCTCCATAAAATGAAAATATCCCGCCCAAGCCGAC
 GTGGGCAAATACACCATGATCTTCCACACCAGAGACAAAGGCAACGAGGTGAACATGGGGCGAATGAAGC
 TGTGCAACAGATGTCCCGGTCTGAAAAACAGATGGCTTATCGAGTTGTTCTTACAGATTACTCTCTGT
 GGAGCACAACCTTTATATGCCAACATCACAGTGGATTCTGGACTGCTGCGTGA

ACGGTACGGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001285793



[View online »](#)

Insert Size:	1035 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001285793.1 , NP_001272722.1
RefSeq Size:	2381 bp
RefSeq ORF:	1035 bp
Locus ID:	56375
UniProt ID:	Q9JJ04
Cytogenetics:	16 B4
Gene Summary:	<p>Responsible for the synthesis of complex-type N-linked oligosaccharides in many glycoproteins as well as the carbohydrate moieties of glycolipids.[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (1) represents the longer transcript. Both variants 1 and 2 encode the same protein. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>